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规格: 50mL/100mL

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**Storage:** 2-8℃ , valid for 1 year.

### Introduction

Trypan Blue Stain Solution, 0.4% is an azo dye, derived from toluidine. It is also known as diamine blue and Niagara blue. It is widely used as a vital stain to distinguish the viable cells from the non-viable cells. It selectively stains the dead tissues or cells blue, while live cells or tissues with intact cell membrane are not colored and are very selective in the compounds that pass through the membrane. So the solution is not absorbed by viable cells. However, it can pass through the cell membrane in a dead cell and so dead cells appear blue under the microscope. Live cells are excluded from staining so this staining method is also described as a dye exclusion method.

Trypan Blue Stain Solution, 0.4% is most commonly used in microscopy for cell counting.

### Protocol(for reference only)

#### 1. Collecting Cells

Centrifuge an aliquot of cell suspension being tested for viability 1 min at 1000-2000 rpm and discard supernatant. Resuspend the cell pellet in PBS.

#### 2. Trypan Blue Stain

Mix 1 part of Trypan Blue Stain Solution, 0.4% and 9 part of cell suspension (cell dilution). Incubate the mixture approximately for 3 mins at room temperature.

#### 3. Cell Count

Take a drop of the trypan blue / cell mixture to a hemacytometer. Place the hemacytometer on the stage of a binocular microscope and focus on the cells. Non-viable cells will be stained blue. Count the unstained (viable) and stained (nonviable) cells separately in the hemacytometer.

$$\frac{\text{total number of viable cells per mL of aliquot}}{\text{total number of cells per mL of aliquot}} \times 100 = \text{viable cells (\%)} \times 100$$

### Note

1. The staining time should not be too long, otherwise the living cells will gradually accumulate dye and stain, which will make the detection result low.
2. If there is precipitation in the dyeing solution before dyeing, it should be filtered out and reused.
3. There is a potential cancer risk. For your safety and health, please wear laboratory clothes and disposable gloves.