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Cat:G1021

Size:2×5mL/2×50mL/2×100mL/2×500mL

Storage: RT, avoid light, valid for at least 2 years.

Kit Components

Reagent	2×5mL	2×50mL	2×100mL	2×500mL	Storage
Reagent(A): WG Stain Solution	5mL	50mL	100mL	500mL	RT,avoid light
Reagent(B): WG Buffer	5mL	50mL	100mL	500mL	RT

Introduction

Wright-Giemsa Stain Kit is a kind of compound dye solution, which has both advantages of Wright and Giemsa stain. It is mainly used in blood and bone marrow smear dyeing. Due to different chemical properties of various cells and components of cells, the affinity of acid dye and basic dye in Wright-Giemsa Stain is also different. Therefore, after the specimen smear is stained with Wright-Giemsa Stain, the corresponding cells show different colors, so as to identify their morphological characteristics. The staining power of Wright stain is strong, which can better show the basophilic degree of the cytoplasm with clear and pure color, but the staining of the nucleus is too deep and the nuclear structure is poor. Giemsa staining has moderate degree of nuclear staining with clear and bright nuclear structure and color, good recognition of nuclear structure, but poor cytoplasmic staining. Therefore, it has the characteristics of good staining effect, clear contrast and simple operation.

Protocol(for reference only)

This kit can be used as the dyeing of many kinds of tissues or cells. Different tissues and cells can be used in different ways. For specific methods, please refer to the existing literature according to your own needs. This solution takes smear as an example to illustrate for reference only.

1. Take the smear and dry it naturally.
2. Drop 2-3 drops of Wright-Giemsa Stain to cover the whole smear, and dye for 1-2 min.
3. Add 0.01M phosphate buffer solution (ph6.4-6.8) in equal amount, shake the slide gently, mix it well with Wright-Giemsa Stain Solution, and dye for 3-5 min.
4. Wash with water, absorb exceed water and view under the microscope.
5. After staining, the cytoplasm and nucleus are clearly stained. The nucleus is purple red with different shades, the cytoplasm is light red, and the particles in the cytoplasm are clearly distinguished.

Note

1. Make smears in appropriate thickness, and fix the smear after drying absolutely, otherwise the cells will fall off easily during the dyeing process.
2. The added dye solution shall not be too small to avoid dye precipitation due to evaporation. Washing time should not be too long to prevent discoloration.
3. Staining is very sensitive to pH, buffer must be used for diluting the dye solution, and the water for washing should be close to neutral, otherwise it may lead to abnormal staining of cells, difficult to identify the morphology, or even errors.
4. If the dyeing is too light, can re-dye. When re-dyeing, add buffer solution first, and then add dye solution. If the dye is too deep, can wash or soak with water or destain with methanol.

Related Products

- P2100 10×Polylysine
 G1010 Giemsa Stain Solution (Working Solution)
 G1040 Wright Stain Solution
 G1100 Eosin Y Stain Solution, Water Solvent
 G1140 Cole's Hematoxylin Solution (For Conventional Stain)