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**Cat:** G1280

**Size:** 2×50ml/2×100ml

**Storage:** 2-8℃, avoid light, valid for 6 months.

### Kit Components

Reagent	2×50ml	2×100ml	Storage
Reagent A: Oxidant	50ml	100ml	2-8℃, avoid light
Reagent B: Schiff Reagent	50ml	100ml	2-8℃, avoid light

### Introduction

In histology, PAS Stain (Periodic Acid-Schiff stain) is mainly used to detect glycogen or other polysaccharides in tissues. The oxidant is periodic acid, which can oxidize the adjacent diol group of polysaccharide into dialdehyde group. The aldehyde group can react with Schiff Reagent to form red insoluble complex.

### Protocol(for reference only)

1. Dewax to distilled water(for paraffin section).
2. Oxidant treat the section for 5-10 min. Rinse in running tap water for 5 min.
3. Stain with Schiff Reagent for 10-15 min. Rinse in running tap water for 5-10 min.
4. (Optional) Re-dyeing with Mayer Hematoxylin Solution. Differentiate and wash; blue and wash.
5. Conventionally dehydrate and transparent by xylene.
6. Seal with resinene.

### Result

Glycogen, Neutral Mucosubstance	Red
Nucleus	Blue

### Note

1. For your safety and health, please wear experimental clothes and disposable gloves.
2. All reagents shall be stored at low temperature, and taken out half an hour before use to restore room temperature.
3. The treatment temperature of Oxidant should not be higher than 20 ℃. When the room temperature is high, the treatment time can be shortened properly. The dyeing time of Schiff Reagent can be adjusted with the temperature. High room temperature can reduce the dyeing time. Low room temperature in winter can extend to about 20 minutes.
4. When using this reagent for the first time, it is recommended to take 1-2 samples for pre-test.

### Related Product

G1010 Giemsa Stain Solution (Working Solution)

G1040 Wright Stain Solution

G1120 Hematoxylin-Eosin/HE Staining Kit

G1080 Mayer's Hematoxylin Solution (For IHC)

G1260 Oil Red O Isopropanol Solution(0.5%)