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## Alkaline Phosphatase Stain Kit (Modified Gomori Ca-Co Method)

**Cat:** G1481

**Size:** 3×50mL

**Storage:** 2-8°C, avoid light, valid for 6 months.

### Kit Components

| Reagent                             | 3×50mL | Storage            |
|-------------------------------------|--------|--------------------|
| Reagent(A): ALP Incubation Solution | 50mL   | 2-8°C, avoid light |
| Reagent(B): Co Solution             | 50mL   | RT, avoid light    |
| Reagent(C): Vulcanizing Solution    | 2×1mL  | RT, avoid light    |
| Reagent(D): ALP Contrast Solution   | 10mL   | 2-8°C, avoid light |

### Introduction

Alkaline phosphatase (ALP or AKP) is a kind of phosphatase widely distributed in mammalian tissues. The optimum pH for its activity is 9.2-9.8. The phosphatase mainly exists in the active sites of substance exchange (cell membrane), such as brush-like margin of intestinal epithelium and proximal convoluted tubule of kidney, stationary cilia of epididymis epithelium, capillary bile duct membrane of liver, and endothelium of arterioles and capillary arteries. It also exists in endoplasmic reticulum, Golgi complex, ingestive vesicles, lysosomes of intestinal epithelium, neutrophils, cell membranes and ingestive vesicles of smooth muscle.

This kit uses metal precipitation method to display the alkaline phosphatase activity. Alkaline phosphatase hydrolyzes to release phosphoric acid using naturally occurring sodium beta-glycerophosphate as the substrate. Then the phosphoric acid is precipitated into calcium phosphate by calcium ion immediately, replaced with cobalt phosphate by Co Solution again, and finally replaced with black precipitation by Vulcanizing Solution.

### Self Provided Materials

Distilled water, Incubator, Water bath

### Protocol (for reference only)

#### For paraffin section staining

1. Dewax paraffin sections and rehydrate in graded alcohol.
2. Incubate the section in ALP Incubation Solution at 37°C for 2-12 h.
3. Rinse with running water for 2mins and place it in distilled water.
4. Incubate the section in Co Solution at 37°C for 5mins.
5. Rinse with running water for 5mins and place it in distilled water.
6. In the above-mentioned process, dilute Reagent (C) to 50 times with distilled water or deionized water to form the Vulcanizing Working Solution (*It is best to prepare it before use and it is not advisable to prepare in advance*). Incubate the section in Vulcanizing Working Solution for 1-2mins.
7. Rinse with running water for 10mins and place it in distilled water.
8. (*Optional*) Re-dying with Nuclear Fast Red and wash with distilled water.
9. Absolute ethanol dehydration, transparent by xylene, seal with resinene.

#### For frozen section staining

1. Fix the frozen section in acetone-chloroform equal mixture at 4°C for 2-5min .
2. Incubate the section in ALP Incubation Solution at 37°C for 45-75min.
3. Rinse with running water for 2min and place it in distilled water.
4. Incubate the section in Co Solution at 37°C for 5min.
5. Rinse with running water for 5mins and place it in distilled water.
6. In the above-mentioned process, dilute Reagent (C) to 50 times with distilled water or deionized water to form the Vulcanizing Working Solution (*It is best to prepare it before use and it is not advisable to prepare in advance*). Incubate the section in Vulcanizing Working Solution for 1-2min.
7. Rinse with running water for 10min and place it in distilled water.
8. (*Optional*) Re-dying with Nuclear Fast Red and wash with distilled water.
9. Seal with Glycerol gelatin.

### Result

|                         |                                     |
|-------------------------|-------------------------------------|
| Positive Site of Enzyme | Black Cobalt Sulfide Precipitation  |
| Nucleus                 | Red(Re-dying with Nuclear Fast Red) |

**Negative Control***(optional)*

1. Reagent (D) is a substrate-free incubation solution. Take the same section into the Reagent (D) - ALP Contrast Solution instead of ALP Incubation Solution, and follow the steps. The negative control result is negative.
2. *(Alternative)*The section can wash in iodine solution and 5% sodium thiosulfate solution each for 3mins and wash with water fully before incubating in incubation solution. This method can use as a negative control.

**Note**

1. ALP Incubation Solution and ALP Vulcanizing Solution are easy to lose effect. It's better to store them in small fractions and use them as soon as open.
2. Vulcanizing Solution is corrosive and operate with caution.
3. When staining frozen sections, the exposure time of sections at room temperature should be reduced.
4. Samples should be fresh and treated immediately after sampling, otherwise the enzyme activity will be affected.
5. Tissue should be fixed in refrigerator at 4°C no more than 24 h, otherwise the activity of enzyme will be weakened or disappeared.
6. For your safety and health, please wear experimental clothes and disposable gloves.