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规格: 4×100mL /4×500mL

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Cat: G1610

Size: 4×100mL /4×500mL

Storage: RT, avoid light, valid for 6 months.

Kit Components

Reagent	4×100mL	4×500mL	Storage
Reagent(A):Hematoxylin Solution	100mL	500mL	RT, avoid light
Reagent(B): Bluing Solution	100mL	500mL	RT
Reagent(C): Orange G Solution	100mL	500mL	RT, avoid light
Reagent(D): EA50 Solution	100mL	500mL	RT, avoid light

Introduction

Papanicolaou method is widely used in routine cytological staining. Papanicolaou stain is initially used to detect estrogen levels in vaginal epithelium and pathogens such as Candida and trichomonad in the genital tract. The combination of orange G6 and EA36 or EA50 can dye the cytoplasm into bright green, blue and pink. At present, most laboratories use finished dye solution, so we should pay attention to the best conditions of each dye solution after improvement. The final cytoplasmic staining should be transparent, and the nuclear chromatin should be easily distinguished.

At present, the improved pasteurization solution contains many kinds of ions, which has polychromatic dyeing efficiency. After staining, the cytoplasm is bright and transparent, and the structure of nuclear membrane, nucleolus and chromatin is clear. The nuclear staining solution is mainly Harris hematoxylin staining solution, and the cytoplasm staining solution is mainly EA36 staining solution and EA50 staining solution. Papanicolaou Stain is used for exfoliated cells. The nucleus is blue or black, and the cytoplasm of keratinized squamous cells is pink or orange red.

Papanicolaou EA50 Stain Kit use EA50 staining solution for cytoplasmic staining and non-toxic modified hematoxylin staining solution for nuclear staining. It is not only suitable for gynecological cytology smear staining, such as screening cervical cancer and precancerous lesions, but also suitable for non gynecological cell samples such as pleural fluid, ascites, sputum and so on.

Self Provided Materials

Fixative(like AF Fixative) ,Series of alcohol, Microscope, Acid alcohol differentiation solution

Protocol(for reference only)

1. For cell smear, fix with 95% alcohol for 10-15mins.
2. Rinse in 95%, 85%, 75% alcohol separately for 1min. Rinse with tap or distilled water for 1min.
3. Stain with Hematoxylin Solution for 5-10mins.Rinse with tap water for 2mins.
4. Differentiate with 1% acid alcohol differentiation solution for about 4-5s or 0.5% acid alcohol differentiation solution for about 10s.Rinse with tap water for 2mins.
5. Blue with Bluing Solution for 2mins.Rinse with tap water for 2mins.
6. Dehydrate in 75%, 85%, 95% alcohol separately for 2mins.
7. Stain with Orange G Solution for 2mins.Dehydrate in 95% alcohol for 2mins twice.
8. Stain with EA50 Solution for 2mins.Dehydrate in 95% alcohol for 2mins twice.
9. Dehydrate in absolute alcohol for 1min twice.
10. Transparent with xylene and seal with resinene.

Result

Nucleus	Blue Purple or Black
Cytoplasm of non keratinized cells	Light Blue or Light Green
Cytoplasm of keratinocytes	Pink or Orange

Note

1. All dye solutions should be filtered before use and changed frequently.
2. For your safety and health, please wear experimental clothes and disposable gloves.