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Cat: G2145

Size: 500mL

Storage: RT, avoid light, valid for 1 year.

### Introduction

The fixative makes the internal structure of the protein change by coagulating and adding compounds, so that the enzyme is inactivated. The fixative can be divided into aldehyde fixative, mercury fixative, alcohol fixative, oxidant fixative, picric acid salt fixative, etc. Formaldehyde in aldehydes and alcohol in alcohols are more commonly used.

Formalin-Zinc (ZF) Fixative, 10% is mainly composed of aldehydes, zinc sulfate and distilled water, which is suitable for the tissue section in most immunohistochemical staining.

### Protocol(for reference only)

Generally, the fixing time of specimens is 1-4h / mm. The fixing time of large specimens should be extended properly.

### Note

1. Formalin-Alcohol-Acid Fixative has certain irritation and corrosiveness. Please operate carefully in ventilated environment
2. The fixing time changes due to different thickness of tissue samples. The suitable thickness of conventional biopsy tissue is 2-4mm, generally no more than 6mm. Proper material selection for the tissue is beneficial to the permeation of the fixative.
3. The effect of temperature on fixation is obvious. Increasing the temperature can accelerate the fixation, but the temperature should not be too high.
4. After fresh tissue is taken out, it should be fixed in time. If it can not be fixed in time, it should be stored in physiological saline for inspection in time.