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Peroxidase Stain Kit(WG-KI Oxidation Method)

Cat: G2371

Size: 2×10mL/2×20mL

Storage: RT, avoid light, valid for 6 months.

Kit Components

Reagent		2×10mL	2×20mL	Storage
Reagent(A): WG-KI Fixative		10mL	20mL	RT, avoid light
Reagent(B): Oxidized WG-KI Solution	B1: WG Solution	1.9mL	3.8mL	RT, avoid light
	B2: KI Solution	0.75mL	1.5mL	RT, avoid light
	B3: WG-KI Buffer	7.5mL	15mL	RT
Before use, mix B1, B1 and B3 as the ratio of 5:2:20 to form Oxidized WG-KI Solution. It is ready to use.				

Introduction

Peroxidase (POX or MPO) is a kind of oxidoreductase produced by microorganisms or plants, which can catalyze many reactions and catalyze the oxidation of the substrate with hydrogen peroxide as the electron acceptor. It mainly exists in the peroxisome of cells. With iron porphyrin as the auxiliary group, it can catalyze the oxidation of phenols and amines by hydrogen peroxide. It has double effects of eliminating the hydrogen peroxide, phenols and amines toxicity.

Peroxidase Stain Kit(WG-KI Oxidation Method) is a modified Pereira POX staining method, which can be used for peroxidase staining of blood, bone marrow or cell smear. The active part of pox is brownish red to blue black particles, which are located in the cytoplasm.

Protocol(for reference only)

1. After the blood and bone marrow smears are dried, drop WG-KI Fixative to cover the blood membrane and fix for 30s.
2. Wash slightly with water, dry in air or absorb the excess solution with filter paper.
3. Add Oxidized WG-KI Solution and dye for 40-60s.
4. Discard the solution and absorb the excess with filter paper. View under the microscope.

Result

Positive Reaction	Brownish Red to Blue Black particles
Negative Reaction	The cytoplasm is blue
Nucleus	Light Blue or Light Purplish Red

In addition to the early neutrophil negative, the positive reaction increases with the maturation of the well-differentiated cells in the following stages, while the senescent neutrophils decreases, monocyte is weak positive and lymphocyte is negative. Plasma cells and megakaryocytes are negative. Eosinophils and Auer bodies show strong positive reaction.

Judgment of positive reaction intensity

Negative	No particle
Weakly Positive	Small particles with sparse distribution
Positive	Coarse particles with dense distribution
Strongly Positive	The particles are coarse, blue black and full of cytoplasm.

Clinical significance

1. In the late stage of acute myeloid leukemia, there are few and large granules.
2. Acute monocytic leukemia cells are negative or weakly positive, there are small and sparse granules.
3. Monocytic acute leukemia is negative or weakly positive.
4. Acute promyelocytic leukemia is strongly positive, some promyelocytes are positive, and malignant cells are negative.
5. Acute lymphocytic leukemia is negative.

Note

1. Blood or bone marrow smear should be fresh, thin and appropriate, fix in time, otherwise it will affect enzyme activity.
2. POX Incubation Solution is easy to lose efficacy or reduce the positive intensity. It is ready to use and not stored for long time.
3. When there are many cells or the positive reaction is weak, the staining time can be appropriately prolonged, but generally it should not be more than 3 minutes. The staining time of AML-M3 cells should be shortened properly.
4. For every staining, the smear of peripheral blood or bone marrow should be taken as negative control.
5. For your safety and health, please wear experimental clothes and disposable gloves.