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Decalcification Endpoint Detection Kit (Chemical Method)

Cat: G2880

Size: 100T

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent	100T	Storage
Reagent (A): Neutralization Solution	100mL	RT
Reagent (B): Ammonium Oxalate Solution	10mL	RT

Introduction

In the process of tissue sectioning, when some tissues contain bone or calcification, the tissue containing calcium should not be directly embedded in paraffin. This is because the density between calcium and paraffin is different, it is difficult to cut a complete section. It is better to fix the calcium containing tissue before decalcification or conduct both at the same time. Then continue operations such as dehydration, transparency, wax immersion, embedding and slicing. There are many decalcification reagents, including organic acid, inorganic acid, EDTA and electrolytic decalcification.

During decalcification, the degree of decalcification should be detected in a while. Excessive decalcification will increase the degree of tissue damage and affect the staining results.

Protocol(for reference only)

1. Take 5mL decalcification solution in use, add 1mL Neutralization Solution, and mix well.
2. Add 0.1mL Ammonium Oxalate Solution and mix well.
3. Observe whether there is precipitation. If there is precipitation, which indicates that decalcification is effective.

Note

1. This method should not be used in the process of decalcification with ion exchange resin.
2. For your safety and health, please wear experimental clothes and disposable gloves.

Appendix:

Determination of the end point of decalcification (physical method): acupuncture, hand pinching, clamp and other methods are used to stop decalcification when the bone tissue becomes soft or there is no sense of resistance during acupuncture. Physical detection will damage the tissue structure to some extent, and try to avoid excessive force or repeated detection.