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## Autophagy/Cytotoxicity Dual Stain Kit (MDC Method)

**Cat:** G0170

**Size:** 100T

**Storage:** -20°C, avoid light, valid for 1 year.

### Kit Components

Reagent	100T	Storage
Reagent A : MDC Stain	1mL	-20°C, avoid light
Reagent B : 10×Wash buffer	20mL	2-8°C
Reagent C : Collection buffer	10mL	2-8°C

### Introduction

Autophagy is a process that cells are stimulated to phagocytize their own cytoplasm or organelles and eventually degrade the phagocytes in lysosomes. Autophagosome is a circular or elliptic structure encapsulated by bilayer membranes. It contains cytoplasm, longevity proteins, complex of abnormal proteins, damages or redundant organelles, viruses and bacteria.

Dansylcadaverine (MDC) is a fluorescent pigment and an eosinophilic stain. It is usually used to detect the formation of autophages as a specific marker dye. Its detection excitation filter wavelength is 355 nm and the blocking filter wavelength is 512 nm. The Autophagy/Cytotoxicity Dual Stain Kit(MDC Method) is suitable for autophagy staining of cultured cells and can counterstain with EB.

### Self Provided Materials

Fluorescence microscope, Low-speed centrifuge, EB stain

### Protocol(for reference only)

#### (一) MDC staining alone

1. Dilute 10× Wash buffer to 1× with distilled water.
2. 800g centrifugation for 5 minutes. Collect the cell in the tube and wash once with 1×Wash buffer of 300-400  $\mu$ l. Remove the supernatant.
3. Add appropriate amount of 1 ×Wash buffer to suspend cells, count and regulate the cell concentration to  $10^6$ /mL.
4. Incubate for 15-45 minutes at room temperature and avoid light.
5. 800g centrifugation for 5 minutes. Collect the cell in the tube and wash twice with 1×Wash buffer of 300-400  $\mu$ l. Remove the supernatant.
6. Add 100 $\mu$ l Collection buffer to suspend the cell, drop it onto the slide and cover the slide.
7. View under the fluorescence microscope (excitation filter wavelength is 355 nm, blocking filter wavelength is 512 nm), count and photograph.

*Note: For adherent cells, can remove the culture medium directly without digestion and stain after washing with 1 × Wash buffer. Adjust the amount of MDC according to the cell number. The dyeing time can be appropriately prolonged and adjusted according to the dyeing results. After dyeing, wash with 1 × Wash buffer and view without adding Collection buffer.*

#### (二) Counterstain with EB

1. Dilute 10× Wash buffer to 1× with distilled water.
2. 800g centrifugation for 5 minutes. Collect the cell in the tube and wash once with 1×Wash buffer of 300-400  $\mu$ l. Remove the supernatant.
3. Add same equal of 1 ×Wash buffer to suspend cells, count and regulate the cell concentration to  $10^6$ /mL.
4. Take 90  $\mu$ l cell suspension into the new tube, add 10  $\mu$ l MDC Stain and 0.2 $\mu$ M EB stain and gently mix.
5. Drop onto the slide and stain it at room temperature for 15-30 minutes without light. Cover the slide.
6. View under the fluorescence microscope (excitation filter wavelength is 512nm), count and photograph.

### Result

Normal cells	Cells are evenly dyed yellow-green fluorescence.
Apoptotic cells	Collapsed chromatin and punctate nucleus are dyed into dense green particles of varying sizes.

### Note

1. MDC Stain and EB Reagent have certain toxicity. Please pay attention to the protection during operation.
2. This kit is often used with acridine orange to distinguish normal cells, apoptotic cells and necrotic cells.
3. During the operation, pay attention to reducing the exposure time of reagents to strong light.
4. For your safety and health, please wear experimental clothes and disposable gloves.