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Nissl Stain Solution (Toluidine Blue Method)

Cat:G1436

Size:100mL

Storage: RT, avoid light, valid for 6 months.

Introduction

Nissl body is a triangular or elliptical block substance distributed in the cytoplasm of nerve cells. It can be dyed purple-blue by methyl violet, methylene blue, toluidine blue and cresyl violet.

Nissl Stain Solution (Toluidine blue method) has many characteristics: simple operation, stable staining, and wide application. It can be used in paraffin tissue sections of Nissl substance/Nissl bodies, neurons and so on. The existence and disappearance of Nissl is an important indicator of whether nerve cells are damaged. When encephalitis, cerebral ischemia and axonal reaction occur, Nissl will dissolve or even disappear.

Protocol(*for reference only*)

1. Fix fresh tissue in 10% formalin solution and routinely dehydrate and embed.
2. Cut paraffin sections in 5 μ m thickness, then dewax paraffin sections and rehydrate in graded alcohol.
3. Stain with Nissl Stain Solution (Toluidine Blue Method) in a temperature box of 50-60 $^{\circ}$ C for 20-40min.
4. Slightly wash with distilled water.
5. Differentiate by 95% ethanol rapidly.
6. Absolute ethanol dehydration, transparent by xylene, seal with resin.

Result

Nissl body	Purple Blue
Nucleus	Light Blue
Background	No Color or Light Blue

Note

1. Nissl dissolves easily in vitro, so the tissue should be fixed immediately after removal, otherwise it is difficult to stain.
2. Tissue fixation plays a very important role. Ethanol or neutral formalin solution can be used for fixation.
3. This staining solution has a good effect on Nissl staining of paraffin tissue sections.
4. Differentiation of 95% ethanol should be carried out rapidly. It is appropriate to observe the slices clearly by naked eye and the background is light blue or colorless.
5. For your safety and health, please wear experimental clothes and disposable gloves.