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Cat: G1533

Size: 4×50mL

Storage: RT, avoid light, valid for 3 months.

Kit Components

Reagent	4×50mL	Storage
Reagent A: Hematoxylin Staining Solution	50mL	RT, avoid light
Reagent B: Acid Alcohol Differentiation Solution	50mL	RT
Reagent C: Sodium Chloride Solution	50mL	RT
Reagent D: Congo Red Staining Solution	50mL	RT, avoid light
Reagent E: Puchtler Alkalization Solution	1mL	RT

Introduction

Amyloid is a kind of extracellular acidophilic substance with no fixed shape, which can exist in different tissues and organs, resulting in diseases called amyloidosis. Amyloid is mainly composed of proteins, most of which are arranged in reverse β - fold structure. Under the electron microscope, the amyloid materials are arranged as fibrils. In the case materials, there are a large number of non branching filaments, most of which are randomly arranged. The histological methods for the identification of amyloid substances include Violet Staining, Congo Red Staining and polarized light microscopy. In 1922, Bennhold found that Congo red can be used to identify starch like substances in vivo, and applied to tissue sections.

Congo Red Amyloid Stain Kit (Puchtler Method) is mainly composed of Congo Red Staining Solution and Hematoxylin Staining Solution. The dyeing principle is that the amyloid has a greater affinity for Congo Red than other tissue structures, and its hydroxyl group combines with the amino group of Congo Red, so that the amyloid is dyed red.

Protocols(for reference only)

1. Conventionally fix , dehydrate and embed.
2. Cut the section in 4 μ m thick, conventionally dewax to distilled water.
3. Soak the section in Hematoxylin Staining Solution for 5mins.
4. Differentiate by Acid Alcohol Differentiation Solution for 2-5s, then immediately remove to distilled water to stop differentiation. Wash twice with water and control to a proper degree under the microscope.
5. Wash with tap water for 2mins.Blue in Scott Bluing Solution or tap water. Wash with tap water for 2mins.
6. Prepare basic Sodium Chloride Solution as the 100:1 ratio of Sodium Chloride Solution:Puchtler Alkalization Solution. It is ready to use. Soak the section in basic Sodium Chloride Solution for 15-20mins.
7. Prepare basic Congo Red Staining Solution as the 100:1 ratio of Congo Red Staining Solution:Puchtler Alkalization Solution. It is ready to use. Soak the section in basic Congo Red Staining Solution for 20-30 mins.
8. Slightly wash with absolute alcohol for 2-3 times.
9. Transparent by xylene and seal with resinene.

Result

Amyloid	Pink Red or Red
Nucleus	Blue

Note

1. Section dewaxing should be as clean as possible, otherwise it will affect the dyeing effect.
2. The acid differentiation solution shall be kept in a closed state and used up as soon as open.
3. Soak dyeing shall be used as much as possible when dyeing with Congo Red Staining Solution. If use drop staining, shall place in a wet box to prevent the solution from volatilizing.
4. basic Sodium Chloride Solution and basic Congo Red Staining Solution are ready to use.
5. For your safety and health, please wear experimental clothes and disposable gloves.

Related Products

G1534 Congo Red Amyloid Stain Kit(Modified Highman Method)

G1530 Congo Red Amyloid Stain Kit, Bennhold

G1535 Congo Red Amyloid Stain Kit(Highman Method)

G1532 Congo Red Amyloid Stain Kit(Modified Stores Method)