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**Cat:** G1594

**Size:** 100mL

**Storage:** 2-8°C, avoid light, valid for 6 months.

### Introduction

In a narrow sense, connective tissue contains three types of fibers: collagen fibers, reticular fibers, elastic fibers. Elastic fibers are found in the lungs, arteries, veins. It is highly refractive, elasticized usually thinner than collagen fibers. Elastic fibers stain well with Gomori aldehyde fuchsin, orcein, resorcin- fuchsin, and Weigert's elastic stain in histological sections.

Aldehyde-Fuchsin Stain Solution can be used for elastic fiber dyeing and other purposes.

### Protocol (for reference only)

1. Fix in 10% neutral formalin fixative, and routinely dehydrate and embed.
2. Cut the paraffin embedded tissue into 4 $\mu$ m-7 $\mu$ m thin sections and routine dewax to water.
3. Operate according to the specific experimental requirements.

### Result

Elastic Fiber	Purple to dark purple
Mast cell granules and mucus	Purple to dark purple
background	Different degrees of yellow

### Note

1. This method can display elastic fiber, anterior elastic fiber and acid resistant fiber, but it needs to slice slightly thicker, 7 $\mu$ m is the best.
2. The oxidizing solution and acidifying solution should not be mixed in advance, and it is better to use them immediately.
3. When aldehyde fuchsin is dyed, it should be capped to prevent the solution from volatilizing.
4. When the Aldehyde-Fuchsin Stain Solution is stored for a long time, the dyeing power will decrease, and the dyeing time should be increased.
5. When islet beta-cells were stained, the time should controlled at 30 min; when pituitary basophils were stained, the time is controlled at 60 min.
6. Orange G dyeing should be light, otherwise it will cover up the color of elastic fiber.
7. For your safety and health, please wear experimental clothes and disposable gloves.