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Ehrlich's Hematoxylin Stain Solution

Cat: G4530

Size: 100mL/500mL

Storage: RT, avoid light, valid for 1 year.

Introduction

Ehrlich's hematoxylin is a natural oxidation hematoxylin like Delafied's. It usually need 3-6 months to finish the process, and the dye ability of the product can be maintained for a long time. Ehrlich hematoxylin, as a good nuclear staining solution, can also be used for mucin staining, such as mucopolysaccharide of cartilage. Ehrlich Hematoxylin Solution is also recommended for bone and cartilage staining.

Ehrlich's Hematoxylin Stain Solution can dye nuclear chromatin clearly and meticulously, and the dyeing time is a little longer, about 15-20min. It is suitable for the production dyeing in teaching and scientific research, but it is not ideal for the frozen section dyeing with this dye solution.

Self Provided Materials

Hydrochloric Acid Ethanol, Bluing Solution, Series Ethanol, Eosin Solution, 4% PFA

Protocol (for reference only)

A. For Paraffin Section Stain

1. Dewaxing To Water
 - ① Dewax with xylene twice and for 5-10min each time.
 - ② (optional) Absolute ethanol acts twice for 3-5 min each time.
 - ③ 95% ethanol for 3-5 min, 90% ethanol for 3-5 min, 80% ethanol for 3-5 min.
 - ④ wash with tap water or distilled water for 1-3min.
2. Staining
 - ① Ehrlich's Hematoxylin Stain Solution for 15-20 min.
 - ② Wash with tap water or distilled water for 5-10S.
 - ③ (optional) Differentiate with Acid Ethanol for 2-5s, wash with tap water for 20-30 s.
 - ④ (optional) Blue with Bluing Solution for 20-40 s, wash with tap water for 30-60 s.
 - ⑤ Stain with Eosin Solution for 10-20s, rinse with tap water for 1-5 s.
3. Dehydration, Transparency and Sealing
 - ① 80% ethanol for 10-20 s, 90% ethanol for 10-20 s, 95% ethanol acts twice for 1-2 min each time, Absolute ethanol acts twice for 2-3 min each time.
 - ② Transparent with xylene for 3 times, 2-3 min each time, Neutral resin seal.

B. For Cell Stain

1. Fix with 4% paraformaldehyde for 10-20 min.
2. Rinse with tap water twice for 2min each time.
3. Wash with distilled water twice for 2min each time.
4. The dyeing, dewaxing, transparency and sealing steps are the same as those of paraffin section, and the action time should be shortened accordingly.

Result

The Nucleus	Blue
the cytoplasm, muscle fiber ,collagen fiber	Recording to the counterstain

Note

1. Slice dewaxing should be as clean as possible. Series ethanol should be replaced frequently.
2. The differentiation time of hydrochloric acid ethanol should be determined according to the slice thickness, tissue type and the old and the new. In addition, the washing time of tap water after differentiation should be enough to wash out the acid thoroughly.
3. 0.2-1% ammonia solution or Scott's blue solution or 0.1-1% lithium carbonate solution are often used as the blueing solution.
4. For your safety and health, please wear experimental clothes and disposable gloves.