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**Cat:** G0110

**Size:** 2 × 100mL

**Storage:** RT, avoid light, valid for 6 months.

### Kit components

Reagent	2 × 100mL	Storage
Reagent(A): Fehling's Solution A	100mL	RT
Reagent(B): Fehling's Solution B	100mL	RT

### Introduction

Fehling's Solution, was invented by German chemist Hermann von Fehling in 1849. Similar to Benedict's reagent, Fehling's Solution is used to detect the presence of reducing sugar. The principle is that the insoluble reducing sugars (glucose, fructose and maltose) can produce brick red cuprous oxide precipitation under the condition of heating. The color change order is: light blue - Brown - brick red (precipitation).

Fehling's Solution is mainly composed of sodium potassium tartrate and copper sulfate. It is often used to identify soluble reducing sugar and urine sugar. It is only used for scientific research and should not be used for clinical testing.

### Self Provided Materials

Test tube and heating device

### Protocols(for reference only)

1. Preparation of Fehling's Working Solution: before use, mix Reagent(A) and Reagent(B) in equal amount to prepare Fehling's Working Solution, which is best to be used immediately.
2. Add 2ml sample to the test tube.
3. Add 1ml Fehling's Working Solution to the test tube and shake well.
4. Place the above mixture in a water bath at 50-5 °C for 2-3min.
5. Observe whether the color of the mixture in the test tube changes. The color change order shall be light blue - brown - brick red (precipitation).

### Result

Reducing sugars (such as ribose, glucose, fructose, etc.)	Brick red precipitation
Non reducing sugar (sucrose, starch, etc.)	No color change

### Note

1. For those with diabetes, it is best to stop taking oral medicine two or three days before the test.
2. This reagent is only used for scientific research and is not suitable for clinical detection.
3. For your safety and health, please wear experimental clothes and disposable gloves.