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Cat: G0980

Size: 10mL/50mL

Storage: 2-8°C, avoid light, valid for 1 year.

Introduction

Mycoplasma Stain Solution(Dienes Method) is a classical reagent which uses azure methylene blue staining to detect Mycoplasma contamination. It is often used to detect apoptosis of cultured adherent or suspended cells and tissue sections. The detection range of cell content by this reagent is generally between 0.1 and 1×10^6 .

Self Provided Material

PBS or Physiological Saline, Slide and Coverslip, Microscope

Protocol(*for reference only*)

(一) Adherent Cells Stain

1. Soak clean coverslip in 70% ethanol for 5 min or longer, dry in sterile super-clean table or wash with sterile PBS or saline for 3 times, and then wash with cell culture medium for 1 time. Place the cover slide in a 6-well plate or other Petri dish. Make the fusion rate till about 50%-80% after overnight incubation.
2. Add the Mycoplasma Stain Solution to cover the surface and incubate for 1-2 min.
3. Remove the solution, wash twice by PBS or physiological saline and 3 min for each time.
4. View under microscope.

(二) Suspended Cells Stain

1. Low-speed centrifugation absorbs most of the liquid and retains about 50 μ l of liquid.
2. Dry slightly and adhere the cell to the slide to make it not easy to flow with the liquid.
3. Add Mycoplasma Stain Solution to cover the surface and incubate for 1-2 min.
4. Remove the dye solution and wash with PBS or physiological saline twice for 3 min each time.
5. View under the microscope.

Result

Mycoplasma	Blue dots were scattered in or around the cells
------------	---

Note

1. In order to obtain the centrifugation of cell precipitation, if the cell precipitation is insufficient, the centrifugal force can be properly increased or prolonged.
2. For your safety and health, please wear experimental clothes and disposable gloves.