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规格: 10mL/100mL/500mL

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Cat: G1066

Size: 10mL/100mL/500mL

Storage: 2-8°C, avoid light, valid for 1 year.

Introduction

The Gram Stain is a different staining technique most widely applied in microbiology. Gram staining is based on the ability of bacteria cell wall to retaining the crystal violet dye during solvent treatment. The cell walls for Gram-positive microorganisms have a higher peptidoglycan and lower lipid content than gram-negative bacteria. Bacteria cell walls are stained by the crystal violet. Iodine is subsequently added as a mordant to form the crystal violet-iodine complex so that the dye cannot be removed easily. However, subsequent treatment with a decolorizer dissolves the lipid layer from the gram-negative cells. As a result, the diffusion of the violet-iodine complex is blocked, and the bacteria remain stained. The length of the decolorization is critical in differentiating the gram-positive bacteria from the gram-negative bacteria. Finally, a counterstain of safranin is applied to the smear to give decolorized gram-negative bacteria a pink color.

Lugol Solution(Gram Method) is a reagent recommended for use in qualitative procedures to differentiate gram-negative from gram-positive organisms.

Protocol(for reference only)

1. Prepare a Slide Smear:

Prepare a thin smear on clear, dry glass slide. Allow to air dry and fix it over a gentle flame, while moving the slide in a circular fashion to avoid localized overheating.

2. Gram Staining:

- 1) Flood with Gram's Crystal Violet Solution for 1 min.
- 2) Wash with tap water.
- 3) Flood the smear with Gram's Iodine Solution for 1 min.
- 4) Pour off the Iodine Solution and gently wash with tap water. Shake off the excess water from the surface.
- 5) Decolorize with Gram's Decolorizing Solution for 20 to 60 s until the blue dye no longer flows from the smear. Further delay will cause excess decolorization in the gram-positive cells, and the purpose of staining will be defeated. Wash with tap water.
- 6) Counterstain with Gram's Safranin Solution for 1 min. Wash with tap water.
- 7) Allow the slide to air dry or blot dry between sheets of clean bibulous paper and view under oil immersion lens.

Result

Gram-positive Organisms	Bluish Purple
Gram-negative Organisms	Pinkish Red

Note

1. The smear of the specimen shall not be too thick and shall be carried out strictly with the operation requirements. If the smear is thick, the decolorization time should be extended until purple no longer appears.
2. The temperature of the slide through the flame shall not be too high.
3. If the Iodine Solution becomes transparent, it cannot be used.
4. When washing with water, the action should be gentle to avoid washing off the bacteria.