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Mayer's Hematoxylin Stain Solution, For IHC

Cat: G1080

Size: 10mL/100mL/500mL

Storage: RT, valid for at least 1 year.

Introduction

Hematoxylin, as an alkaline dye and one of the two main staining reagents for HE staining, is usually used for nuclear staining. It can be roughly divided into two types according to different dyeing operations and applicable samples. One is the progressive dyeing hematoxylin with less hematoxylin content, which is suitable for lining dyeing and usually does not need differentiation. The other is the degenerative staining hematoxylin with high content of hematoxylin and clear and long-lasting staining but needs differentiation.

This product is a progressive hematoxylin staining solution, which is suitable for counterstaining of tissue sections in immunohistochemistry. Usually, it does not need differentiation, and the nuclei are blue after staining. This product is working fluid and can be used directly. The dyeing solution can be reused for many times, but the dyeing effect will gradually decrease.

Protocol(for reference only)

1. Rinse sections in distilled water to remove chromogen/substrate solution from immunohistochemically stained tissue or cell.
2. Incubate in Mayer's Hematoxylin Stain Solution (For IHC) for 30s-2min or 5-10min depending on desired signal intensity.
3. (optional) Differentiate with differentiation solution for 30s.
4. Wash with tap water for 10min or other anti blue solution for 3-10min to fully anti blue.
5. Dehydrate in gradient alcohol, transparent by xylene and seal with resinene.

Result

Nucleus	Blue
---------	------

Note

1. When using this kit for the first time, it is recommended to take 1-2 samples for pre-test.
2. The best temperature is 2-8 °C, which is at least one year.
3. Light staining is recommended in the process of staining. Generally, it is only necessary to be able to distinguish the nucleus. Too deep color may affect the color of the cytoplasm.
4. For your safety and health, please wear experimental clothes and disposable gloves.