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Cat:G1261

Size:2×50mL/2×100mL

Storage:2-8°C, avoid light, valid for 6 months.

Kit Components

Reagent		2×50mL	2×100mL	Storage
Reagent(A):Modified Oil Red O Stain Solution	A1: Oil Red O Stain A	30mL	60mL	2-8°C,avoid light
	A2: Oil Red O Stain B	20mL	40mL	RT
Mix A1 and A2 with the ratio of 3:2 and place for 10 min to form Modified Oil Red O Stain Solution. Configure before use.				
Reagent(B):Mayer's Hematoxylin Solution		50mL	100mL	2-8°C,avoid light

Introduction

Neutral fat stains often use Sudan II, Sudan III, Sudan IV, Sudan black B, oil red O and so on. Sudan stain is often used in traditional methods. Recently, azo dye oil red O is more suitable for dyeing fat. Oil red O is a strong lipid solvent and dye. It is easy to bind to triglycerides in droplet shape, but has a slightly poor binding to phospholipids. The dyeing principle is generally regarded as a physical miscibility or adsorption, and the fat is dyed by miscibility. The solubility of dyes in frozen sections is higher than that in the original solvent, so when dyeing, the dyes are transferred from organic solvents to lipids and the fats are dyed.

Modified Oil Red O Stain Solution is mainly used to show fatty degeneration of tissues and organs and abnormal lipid-like sedation. It often occurs in fatty degeneration of liver, kidney, heart and other parenchymal organs that most of the neutral fat droplets appear in cells. It can identify and diagnose tumors and their properties in adipose tissue. The sample should not be fixed with ethanol (10% formalin if required) or paraffin. Frozen sections or carbon wax sections were needed. The positive staining results of fat were orange to red, but the specific color depended on the concentration of lipid.

Self Provided Materials

60% isopropanol, distilled water, 1% hydrochloric acid solution, dilute lithium carbonate solution, glycerol gelatin or arabinose

Protocol(for reference only)

1. Cut frozen sections 6-10 um thick,unfix or wash by water after fix in 10% formalin for 10min.
2. Rinse sections slightly with distilled water.
3. Soak sections in 60% isopropanol for 20-30 s.
4. Stain section in Modified Oil Red O Stain Solution (capped) for 10-15 min.
5. Colour separation: Wash slightly in 60% isopropanol to remove the dye solution.
6. Rinse sections slightly with distilled water.
7. Re-dyeing by Mayer's Hematoxylin Solution for 1-2 min.
8. (optional)Slightly differentiate by 1% hydrochloric acid solution.
9. (optional)Rinse by tap water for 10 min or blue in dilute lithium carbonate solution.
10. Rinse sections slightly with distilled water.

Result

Neutral Fat	Orange or Red
Nucleus	Blue

Note

1. The Modified Oil Red O Stain Solution is not stable enough and easy to precipitate, so it is not suitable to prepare it in advance.
2. If 60% isopropanol is not available, 70% ethanol can also be used.
3. Because fat is soluble in organic solvents, it is shown that fat can not be treated like paraffin slices, but can be shown by frozen section staining.
4. The frozen sections stained with fat should not be too thin, and too thin sections often cause lipid loss.
5. Mayer's Hematoxylin Solution should not stain for too long.

6. The dyeing result can not last for a long time, and should be viewed and photographed as soon as possible.
7. Samples sealed with glycerol gelatin can not last for a long time. If long-term preservation is required, it is recommended to seal with resinene between the cover slide and the slide.
8. For your safety and health, please wear experimental clothes and disposable gloves.