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Storage: 2-8°C, valid for 6 months.

Introduction

Glycogen staining is one of the conventional staining methods in pathology. McManus first used PAS technology to display mucin in 1946. This method is often used to display glycogen and other polysaccharides. The staining solution can not only display glycogen, but also show neutral mucilaginous substances and some acidic substances, as well as cartilage, pituitary, mould, fungus, pigment, amyloid substance, basement membrane, etc.

Diastase Aqueous Solution, 1%, pH5.3 is mainly used for section treatment before glycogen PAS staining. When glycogen is digested, prepare two identical sections. After dewaxing, one section is treated with Diastase Solution, 1%, pH5.3, Water Solvent, and the other one is treated with PBS or distilled water. Then the two sections are stained with PAS method, and the disappearance of staining after digestion indicated the presence of glycogen.

Protocol(for reference only)

1. Prepare two identical sections and dewax to distilled water.
2. Treat one section into Diastase Solution, 1%, pH5.3, Water Solvent at 37°C for 20min. The other one directly into distilled water for 1h without treated by Diastase Solution, 1%, pH5.3, Water Solvent.
3. Rinse the two sections in tap water respectively for 5-10min.
4. Follow the glycogen PAS staining steps.

Result

Glycogen, Neutral and Salivary Mucin	Red to Purple
Various Glycoprotein	Red to Purple
Nucleus	Blue
The section untreated , glycogen is bright red or red purple; The section treated by Diastase Solution, glycogen shows negative.	

Note

1. Section dewaxing should be as clean as possible, otherwise it will affect the dyeing effect.
2. A positive control section is required to verify the enzyme activity.
3. Avoid too much sunlight and air during use. Before use, it is better to take it out in advance and restore room temperature, and use it in dark.
4. The staining time of frozen section should be as short as possible.
5. For your safety and health, please wear experimental clothes and disposable gloves.