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Neutral Red Sterile Solution,0.33%, For Cells

Cat: G1310

Size: 100mL/500mL

Storage: 2-8°C, valid for 1 year.

Introduction

Neutral Red Sterile Solution,0.33%, For Cells is a kind of common staining solution which can dye the nucleus red for tissue or cell staining. This staining solution is prepared with neutral pH isotonic buffer and can be directly used for staining of living cells or tissues.

The staining solution also can be directly used for the staining of cultured cells. The nucleic acid in the nucleus is acidic, so the nucleus will be dyed red. Because lysosomes are also in acidic environment, lysosomes can also be dyed red by Neutral Red Solution.

Protocol(*for reference only*)

1. Sample Treatment

Remove the medium and wash once with appropriate buffer such as PBS, DPBS or Hanks.

2. Neutral Red Solution Staining

Remove the washing solution, add a proper amount of staining solution to ensure that the cells are fully covered and stained for 2-10 min (the time can be adjusted according to the staining results and requirements). Wash with PBS, DPBS or Hanks for 1-2 times, and then observe and photograph.

Note

1. When using this reagent for the first time, it is recommended to take 1-2 samples for pretest.
2. For your safety and health, please wear experimental clothes and disposable gloves.
3. Long term storage of Neutral Red Solution will produce precipitation. It can be directly used to absorb the supernatant of the dye solution, or it can be used to remove the precipitate by filtration or centrifugation, without affecting the use effect.