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Cat: G1316

Size: 100mL/500mL

Storage: 2-8°C, avoid light, valid for 1 year.

Introduction

Neutral red acts as a pH indicator, changing from red to yellow between pH 6.8 and 8.0. The molecular formula is $C_{15}H_{16}N_4 \cdot HCl$, and the molecular weight is 288.78. In neutral or slightly alkaline environment, living cells of plants can absorb a large number of neutral red and excrete into vacuoles. Neutral Red Solution, 0.1% can be used as a red nuclear counterstain in various histological procedures.

Protocol (for reference only)

1. Sample Treatment

For paraffin section, dewax in xylene twice for 5-10mins each. Rehydrate with series ethanol (100%, 90%, 80%, 70%) for each level 3mins, and finally in distilled water 3min.

For frozen section restore to RT and wash with distilled water for 2mins.

For cultured cell, fix with 4% PFA for more than 10mins. Then wash with distilled water twice for 2mins each.

2. Neutral Red Solution Staining

- a. For the above treated samples, stain with Neutral Red Solution, 0.1% for 0.5-1h.
- b. Rinse briefly in distilled water or distilled water.
- c. View under the microscope and photograph.

Note

1. If dehydration, transparency and sealing treatment are needed, xylene, resinene or other sealing agent shall be prepared by oneself. If the sample is paraffin section, anhydrous ethanol and xylene shall be prepared by ourselves. The color may fade during dehydration. It is recommended to observe visually.
2. When using this reagent for the first time, it is recommended to take 1-2 samples for pre-test.
3. When staining with Neutral Red Solution, the time can be adjusted according to the dyeing results and requirements.
4. For your safety and health, please wear experimental clothes and disposable gloves.