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Cat: G1340

Size: 7×50mL/7×100mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent	7×50mL	7×100mL	Storage	
Reagent (A): Weigert's Iron Hematoxylin Solution	A1: Weigert's Solution A	25mL	50mL	RT, avoid light
	A2: Weigert's Solution B	25mL	50mL	RT, avoid light
Mix equal parts of A1 and A2 to form Weigert's Iron Hematoxylin Solution, which is stable for about 4 h.				
Reagent (B): Acid Alcohol Differentiation Solution	50mL	100mL	RT	
Reagent (C): Bluing Solution	50mL	100mL	RT	
Reagent (D): Ponceau-Acid Fuchsin Solution	50mL	100mL	RT, avoid light	
Reagent (E): Acetic Acid Solution	50mL	100mL	RT	
Reagent (F): Phosphomolybdic Acid Solution	50mL	100mL	RT, avoid light	
Reagent (G): Aniline Blue Solution	50mL	100mL	RT, avoid light	

Introduction

Masson Trichrome Stain Kit is intended for use in the study of connective tissue, muscle and collagen fibers. It is mainly used in distinguishing collagen from smooth muscle since these two components look similar under the microscope. In the Masson's Trichrome Stain Kit nuclei are stained with Weigert's iron hematoxylin, and cytoplasm and muscle are then stained with Beibrich scarlet-acid fuchsin. After treatment with phosphotungstic and phosphomolybdic acid, collagen is demonstrated by staining with aniline blue. Rinsing in acetic acid after staining renders the shades of color more delicate and transparent.

This kit may be used on formalin-fixed, paraffin-embedded or frozen sections.

Protocol (for reference only)

- Dewax to distilled water.
- Stain with Weigert's Iron Hematoxylin Solution for 5-10mins.
- Differentiate with Acid Alcohol Differentiation Solution for 10-15s.
- Blue in Bluing Solution for 2-5mins. Rinse in deionized water.
- Stain with Ponceau-Acid Fuchsin Solution for 5-10mins. Rinse in deionized water.
- Differentiate in Phosphomolybdic Acid Solution for 1-2mins or till collagen is not red.
- Without rinsing, add Aniline Blue Solution to section and stain for 1-2mins.
- Place section in Acetic Acid Working Solution (Mix 1 part of Acetic Acid solution and 2 part deionized water) for 1 min.
- Dehydrate very quickly in 95% ethanol, absolute ethanol (these steps will wipe off Beibrich Scarlet-Acid Fuchsin Staining) and transparent in xylene.
- Seal with resinene.

Result

Nucleus	Black
Cytoplasm, Muscle fibers	Red
Collagen	Blue

Note

- Slice dewaxing should be as clean as possible.
- Equal amount of A1 and A2 were mixed to form Weigert iron hematoxylin staining solution, which generally lost its dyeing ability within 24 h.
- Tissue fixation plays a very important role. Different fixatives can prolong or shorten the dyeing time.
- In the classic Masson trichromatic staining, Harris hematoxylin was used to stain the nucleus, but the color of the section after Harris hematoxylin staining was not bright enough. Weigert hematoxylin was used to stain the nucleus in this staining solution, because the purpose of staining is mainly to distinguish collagen fibers and muscle fibers, and this staining step can also be omitted generally.

5. The differentiation time of acid ethanol should be determined according to the thickness of slice, the type of tissue and the old and new.
6. Weak acid solution can make the color clearer and brighter. If it is used in a large amount, it can be replaced by 0.1-0.3% acetic acid solution. If aniline blue is over dyed, the weak acid cleaning steps in steps 6 and 7 can be omitted, and aniline blue dyeing can be carried out directly after phosphomolybdic acid cleaning.
7. The differentiation of phosphomolybdic acid should be controlled under microscope until the collagen fibers are light red and the fibers are red. The differentiation time depends on the depth of staining, generally 1-2mins.
8. Masson Bluing Solution can also be replaced by Scott bluing solution or 0.1-1% lithium carbonate solution.