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Masson's Trichrome Stain Kit (Fast Green FCF Method)

Cat: G1343

Size: 7×50mL/7×100mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent		7×50mL	7×100mL	Storage
Reagent(A) :Weigert Hematoxylin Dyeing Solution	A1:Weigert Stain Solution A	25mL	50mL	RT, avoid light
	A2:Weigert Stain Solution B	25mL	50mL	RT, avoid light
Before use, mix with A1 and A2 in equal amount to form the Weigert Hematoxylin Dyeing Solution. It is not suitable to prepare it in advance.				
Reagent(B): Acid Ethanol Differentiation Solution		50mL	100mL	RT
Reagent(C): Masson Bluing Solution		50mL	100mL	RT
Reagent(D): Ponceau Fuchsin Solution		50mL	100mL	RT, avoid light
Reagent(E): Weak Acid Solution		50mL	100mL	RT
Reagent(F): Phosphomolybdic Acid Solution		50mL	100mL	RT, avoid light
Reagent(G): Fast Green Solution		50mL	100mL	RT, avoid light

Introduction

In a narrow sense, connective tissue contains three kinds of fibers: collagen fibers, reticular fibers and elastic fibers. Collagen fibers are the most widely distributed and abundant fibers. Masson's Trichrome Staining is the most classical method in connective tissue staining. It is an authoritative and classical technique for collagen fibers staining. The Masson's Trichrome Staining usually refers to the staining of nucleus and the selective display of collagen and muscle fibers. The dyeing principle of this kit is related to the size of anionic dye molecules and the permeability of the tissues: the size of the molecules is reflected by the molecular weight, and the small molecular weight can easily penetrate the compact and low permeability tissues, while the large molecular weight can only enter the loose and high permeability tissues. However, the molecular weight of light green or aniline blue is very large, so after Masson's Trichrome Staining, muscle fibers are red, collagen fibers are green (light green) or blue (aniline blue), mainly used to distinguish collagen fibers from muscle fibers.

The characteristics of Masson's Trichrome Staining are: (1) stable dyeing; (2) short differentiation time, 1-2 s; (3) clear and bright color; (4) wide application, suitable for tissue paraffin section, frozen section and other dyeing; (5) long preservation time and not easy to fade.

Self Provided Materias

Fixative:10% formalin or 10% Neutral formalin., Series of ethanol, Distilled water.

Protocol(for reference only)

1. Dewax paraffin sections and rehydrate in graded alcohol.
2. Stain with Weigert Hematoxylin Dyeing Solution for 5-10mins.
3. Differentiat by Acid Ethanol Differentiation Solution for 1-2 seconds, then wash with water.
4. Bluing by Masson Bluing Solution for 3-5mins, then wash with water.Wash with distilled water for 1min.
5. Stain with Ponceau Fuchsin Solution for 5-10mins.
6. In the above operation process, mix distilled water and Weak Acid Solution with the ratio of 2:1 to form Weak Acid Working Solution . Wash with Weak Acid Working Solution for 1min.
7. Wash with Phosphomolybdic Acid Solution for 1-2mins.Wash with Weak Acid Working Solution for 1min.
8. Stain with Fast Green Solution for 1-2mins.Wash with Weak Acid Working Solution for 1min.
9. Quick dehydration of 95% ethanol. Dehydrated by Anhydrous ethanol for three times, each time for 5 to 10 s.
10. Transparentize by xylene for three times and each time for 1-2mins. Seal with resinene.

Result

Collagen Fiber/Protein	Green
Cytoplasm, Muscle And Red Blood Cell	Red
Nucleus	Bluish Brown

Note

1. Slice dewaxing should be as clean as possible. Fixation plays an important role, and dyeing time can be prolonged or shortened by using different fixatives.
2. Weigert Hematoxylin Dyeing Solution is obtained by mixing A1 and A2 equally. It is generally recommended to be used within 4 hours.
3. The differentiation time of Acid Ethanol Differentiation Solution should be determined by the thickness of sections, the type of tissues and the old and new ones.
4. Weak Acid Solution can make the color clearer and brighter. If use in large quantities, 0.1-0.3% acetic acid solution can be prepared to replace it.
5. The differentiation of Phosphomolybdic Acid Solution should be controlled under the microscope until that the collagen fibers are light red and the fibers are red. Differentiation time depends on the depth of staining, generally for 1-2mins.
6. Masson Bluing Solution can also be replaced by Scott Blue Promoting Solution or 0.1-1% Lithium Carbonate Aqueous Solution.