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- [1] Dan Xiao, Yue Zhang, Rui Wang, et al. Emodin alleviates cardiac fibrosis by suppressing activation of cardiac fibroblasts via upregulating metastasis associated protein 3. *Acta Pharmaceutica Sinica B*. April 2019. (IF 5.808)
- [2] Na Li, Lin Zhou, Weilong Xie, et al. Alkaline phosphatase enzyme-induced biomineralization of chitosan scaffolds with enhanced osteogenesis for bone tissue engineering. *Chemical Engineering Journal*. September 2019;371:618-630. (IF 8.355)
- [3] Ke Xue, Jun Zhang, Cong Li, et al. The role and mechanism of transforming growth factor beta 3 in human myocardial infarction-induced myocardial fibrosis. *Journal of Cellular and Molecular Medicine*. April 2019. (IF 4.658)

Modified Masson's Trichrome Stain Kit

Cat: G1346

Size: 8×50mL/8×100mL

Storage: 2-8°C, avoid light, valid for 1 year.

Kit Components

Reagent	8×50mL	8×100mL	Storage
Reagent(A): Mordant Solution	50mL	100mL	RT
Reagent(B): Celestite Blue Solution	50mL	100mL	2-8°C, avoid light
Reagent(C): Mayer Hematoxylin Solution	50mL	100mL	2-8°C, avoid light
Reagent(D): Acid Differentiation Solution	50mL	100mL	RT
Reagent(E): Ponceau-Acid Fuchsin Solution	50mL	100mL	RT, avoid light
Reagent(F): Phosphmolybic Acid Solution	50mL	100mL	RT, avoid light
Reagent(G): Aniline Blue Solution	50mL	100mL	RT, avoid light
Reagent(H): Acetic Acid Solution	50mL	100mL	RT

Introduction

Masson Trichrome Stain is intended for use in the study of connective tissue, muscle and collagen fibers. It is mainly used in distinguishing collagen from smooth muscle since these two components look similar under the microscope.

The difference between the Modified Masson's Trichrome Stain Kit and the conventional Masson Trichrome Stain Kit is that the nucleus of the cell is light stained with Celestite Blue Hematoxylin Solution. It has the following characteristics: short differentiation time, clear and bright colors; wide application range, suitable for paraffin section, frozen section and other staining of tissues; long preservation time of stained sections and not easy to fade. Modified Masson's Trichrome Stain Kit shows that collagen fiber is blue, muscle fiber, cytoplasm, cellulose, keratin and erythrocytes are red, and nucleus is blue, which was mainly used to distinguish collagen fiber and muscle fiber.

Self Provided Materials

10% Formalin, Distilled water, Series of ethanol, xylene, Dye vat.

Protocol (for reference only)

1. Fix tissues in 10% formalin, then dehydrate and embed.
2. Cut the section in 4µm thick, dewax to distilled water.
3. Incubate the section into Mordant Solution in 37°C overnight or incubate in 57-60°C incubator for 1h, then wash with running water for 10mins. (See note 2)
4. Drop Celestite Blue Solution onto the section and stain for 2-3mins. Slightly wash with distilled water twice, each for 15-30s.
5. Drop Mayer Hematoxylin Solution onto the section and stain for 2-3mins. Slightly wash with distilled water twice, each for 15-30s.
6. Differentiate by Acid Differentiation Solution for several seconds. Rinse in running water for 10mins. (See note 2)

7. Drop Ponceau-Acid Fuchsin Solution onto the section and stain for 10mins. Slightly wash with distilled water twice,each for 15-30s .
8. Treat with Phosphmolybic Acid Solution for 10mins.
9. Discard the remaining dye solution and directly stain with Aniline Blue Solution for 5mins without washing with water.
10. As the radio of 1:2, mix the Acetic Acid Solution and water to prepare Acetic Acid working Solution, and rinse in Acetic Acid working Solution for 2mins.
11. Dehydrate quickly in 95% ethanol. Dehydrate in absolute ethanol 3 times for each time 5-10s.
12. Transparent by xylene three times for each time 1-2mins. Seal with resinene.

Result

Collagen Fiber	Blue
Muscle Fiber, Cytoplasm, Cellulose, Keratin and Erythrocytes	Red
Nucleus	Bluish Brown

Note

1. Section dewaxing should be as clean as possible.
2. The Mordant Solution and Acid Differentiation Solution of this product are optimized formulas, which are colorless solvents without dangerous chemicals. They are quite different from the properties of classical picric acid mordant and acid ethanol differentiation solution, and belong to normal phenomena. The above two reagents are optimized substitutes for the original reagent. If you prefer the classic formula, you can replace the reagents used in the corresponding steps by yourself.
3. The differentiation time of Acid Differentiation Solution should be determined by the thickness of section, the type of tissue and the old and new.
4. Phosphmolybic Acid Solution has two effects. One is that Phosphmolybic Acid Solution can make the collagen fiber dyed red differentiate into colorless or light red, while the muscle fiber cellulose is still bright red; the other one is that it can also play a role of mordant dyeing on the collagen fiber, so that the collagen fiber is easy to combine with Aniline Blue Solution of macromolecular dye.
5. Treat the section with Acetic Acid Solution after Aniline Blue Solution staining to remove the blue in the original pulp and make the dyeing bright and clear. If the tissue is fixed with Zenker solution, the treatment with Acetic Acid Solution can be prolonged to 5 mins.
6. Acetic Acid Solution can make the color more clear and bright. If the use amount is large, can replace it with 0.1-0.3% acetic acid solution by yourself .
7. For your safety and health, please wear experimental clothes and disposable gloves.

Reference

- [1] Dan Xiao,Yue Zhang,Rui Wang,et al. Emodin alleviates cardiac fibrosis by suppressing activation of cardiac fibroblasts via upregulating metastasis associated protein 3. *Acta Pharmaceutica Sinica B*. April 2019. (IF 5.808)
- [2] Na Li,Lin Zhou,Weilong Xie,et al. Alkaline phosphatase enzyme-induced biomineralization of chitosan scaffolds with enhanced osteogenesis for bone tissue engineering. *Chemical Engineering Journal*. September 2019;371:618-630. (IF 8.355)
- [3] Ke Xue,Jun Zhang,Cong Li,et al. The role and mechanism of transforming growth factor beta 3 in human myocardial infarction-induced myocardial fibrosis. *Journal of Cellular and Molecular Medicine*. April 2019. (IF 4.658)