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Nissl Stain Kit (Methyl Violet Method)

Cat:G1432

Size:2×50mL/2×100mL

Storage:RT, avoid light, valid for 6 months.

Kit Components

Reagent	2×50mL	2×100mL	Storage
Reagent(A):Methyl Violet Solution	50mL	100mL	RT, avoid light
Reagent(B): Nissl Differentiation	50mL	100mL	RT

Introduction

The neuronal soma consists of a large nucleus with a folded nuclear membrane, a sparse chromatin and a distinct nucleolus. In the cell body, the cytoplasm is Nissl granules, which can represent rough endoplasmic reticulum and produce specific punctate basophilic granules in many neurons. Nissl granules can be shown by many dyes, such as neutral red, methylene blue, toluidine blue and methyl violet. Variation of dyeing, pH and differentiation time make some staining not only highlight Nissl substance, but also include the nucleus and glia cells. Nissl body is a triangular or elliptical block substance distributed in the cytoplasm of nerve cells. It can be dyed purple-blue by methyl violet, methylene blue, toluidine blue and tar violet.

This kit has many characteristics: simple operation, stable staining, short differentiation time and wide application. It can be used in paraffin tissue sections of Nissl substance/Nissl bodies, neurons and so on.

Protocol (for reference only)

1. Fixation: fix in ethanol, Carnoy fixative or neutral formalin salt solution.
2. Tissue section: the thickness of paraffin section is 7-10 or 25 μ m (see note 4).
3. Dewax paraffin sections and rehydrate in graded alcohol.
4. Stain with Methyl Violet Solution for 10-20min.
5. Rinse with distilled water.
6. Differentiate with Nissl Differentiation for 4-8 s until most of the staining solution is eliminated.
7. Absolute alcohol dehydration, transparentize by xylene, view under the microscope.
8. If necessary, repeat step 6 and 7, and when repeating, add a small amount of Nissl Differentiation .
9. Rinse thoroughly in xylene.
10. Seal with Canadian balm or DPX.

Result

Nissl body	Violet Blue to Violet Black
Neuron	Light Violet Blue
Nucleus	Violet Blue

Note

1. Nissl dissolves easily in vitro, so the tissue should be fixed immediately after removal, otherwise it is difficult to stain. Tissue fixation plays a very important role. Ethanol, Carnoy or neutral formalin solution can be used for fixation.
2. This staining solution has a good effect on Nissl staining of paraffin tissue sections.
3. If only to confine the Nissl substance, it must be differentiated by acid differentiation solution after staining.
4. Paraffin section thickness is 7-10 μ m or 25 μ m (cortical neuron density should be assessed with 25 μ m section).
5. The stained specimens must be kept away from light, otherwise they will fade easily.
6. For your safety and health, please wear experimental clothes and disposable gloves.