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## Congo Red Amyloid Stain Kit (Bennhold Method)

**Cat:** G1530

**Size:** 5×50mL

**Storage:** RT, avoid light, valid for 1year.

### Kit Components

Reagent	5×50mL	Storage
Reagent A: Bennhold Hematoxylin Staining Solution	50mL	RT, avoid light
Reagent B: Acid Differentiation Solution	50mL	RT
Reagent C: Scott Bluing Solution	50mL	RT
Reagent D: Congo Red Staining Solution	50mL	RT, avoid light
Reagent E: Bennhold Differentiation Solution	50mL	RT

### Introduction

Amyloid is a kind of extracellular acidophilic substance with no fixed shape, which can exist in different tissues and organs, resulting in diseases called amyloidosis. Amyloid is mainly composed of proteins, most of which are arranged in reverse  $\beta$  - fold structure. Under the electron microscope, the amyloid materials are arranged as fibrils. In the case materials, there are a large number of non branching filaments, most of which are randomly arranged. The histological methods for the identification of amyloid substances include Violet Staining, Congo Red Staining and polarized light microscopy. In 1922, Bennhold found that Congo red can be used to identify starch like substances in vivo, and applied to tissue sections.

Congo Red Amyloid Stain Ki (Bennhold Method) is mainly composed of Congo Red Staining Solution and Bennhold Hematoxylin Staining Solution. The dyeing method has stable performance and is a very classical method for dyeing amyloid.

### Protocols (for reference only)

1. Conventionally fix, dehydrate and embed.
2. Cut the section in 4 $\mu$ m thick, conventionally dewax to distilled water.
3. Soak the section in Bennhold Hematoxylin Staining Solution for 5mins.
4. Differentiate by Acid Alcohol Differentiation Solution for 2-5s, then immediately remove to distilled water to stop differentiation. Wash with tap water for 2mins.
5. Blue in Scott Bluing Solution. Wash with tap water for 2mins.
6. Soak the section in Congo Red Staining Solution for 30mins.
7. Quickly differentiate by Bennhold Differentiation Solution. Wash with tap water for 1-2mins.
8. Dehydrate by series of alcohol, transparent by xylene, then seal with resinene.

### Result

Amyloid	Red
Nucleus	Blue

### Note

1. Section dewaxing should be as clean as possible, otherwise it will affect the dyeing effect.
2. The acid differentiation solution shall be kept in a closed state and used up as soon as open.
3. Soak dyeing shall be used as much as possible when dyeing with Congo Red Staining Solution. If use drop staining, shall place in a wet box to prevent the solution from volatilizing.
4. The differentiation process of Bennhold Differentiation Solution is very important. Pay attention to the degree of differentiation until there is no excess dye flowing down. When the differentiation time is short, the collagen fiber is also dyed red; when the differentiation is excessive, the amyloid is also decolorized. If the decolorization is excessive, can clean the section and then soak with Congo Red Staining Solution again.
5. For your safety and health, please wear experimental clothes and disposable gloves.

### Related Products

*G1534 Congo Red Amyloid Stain Kit (Modified Highman Method)*

*G1533 Congo Red Amyloid Stain Kit (Puchtler Method)*

*G1532 Congo Red Amyloid Stain Kit (Modified Stores Method)*