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Congo Red Amyloid Stain Kit(Modified Highman Method)

Cat: G1534

Size: 3×50mL

Storage: RT, avoid light, valid for 6 months.

Kit Components

Reagent	3×50mL	Storage
Reagent A:Modified Highman Solution	50mL	RT, avoid light
Reagent B:Highman Differentiation Solution	50mL	RT
Reagent C:Mayer Hematoxylin Staining Solution	50mL	RT, avoid light

Introduction

Amyloid is a kind of extracellular acidophilic substance with no fixed shape, which can exist in different tissues and organs, resulting in diseases called amyloidosis. Amyloid is mainly composed of proteins, most of which are arranged in reverse β - fold structure. Under the electron microscope, the amyloid materials are arranged as fibrils. In the case materials, there are a large number of non branching filaments, most of which are randomly arranged. The histological methods for the identification of amyloid substances include Violet Staining, Congo Red Staining and polarized light microscopy. In 1922, Bennhold found that Congo red can be used to identify starch like substances in vivo, and applied to tissue sections.

Congo Red Amyloid Stain Kit(Modified Highman Method) also known as Methanol Congo Red Staining, which is mainly composed of Modified Highman Solution and Mayer Hematoxylin Staining Solution.

Protocols(for reference only)

1. Conventionally fix in 10% neutral formalin, dehydrate and embed.
2. Cut the section in 4 μ m thick, conventionally dewax to distilled water.
3. Soak the section in Modified Highman Solution for 10mins and discard the excess solution.
4. Differentiate by Highman Differentiation Solution for 2-5s, then immediately remove to distilled water to stop differentiation. Wash twice with water and control to a proper degree under the microscope.
5. Wash with tap water for 5mins.
6. Dyeing with Mayer Hematoxylin Staining Solution for 1-2mins or less time.Wash with tap water for 5mins.
7. Dehydrate by series of alcohol, transparent by xylene, then seal with resinene.

Result

Amyloid	Red
Nucleus	Blue

Note

1. Section dewaxing should be as clean as possible, otherwise it will affect the dyeing effect.
2. The acid differentiation solution shall be kept in a closed state and used up as soon as open.
3. Soak dyeing shall be used as much as possible when dyeing with Modified Highman Solution. If use drop staining, shall place in a wet box to prevent the solution from volatilizing.
4. The differentiation process of Highman Differentiation Solution is very important. When the differentiation time is short, the collagen fiber is also dyed red; when the differentiation is excessive, the amyloid is also decolorized. If the decolorization is excessive, can clean the section and then soak with Modified Highman Solution again.
5. Dehydration should be rapid to avoid decolorization.
6. For your safety and health, please wear experimental clothes and disposable gloves.