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**Cat:** G1541

**Size:** 2×100mL

**Storage:** RT, avoid light, valid for 1 year.

### Kit Components

Reagent	2×100mL	Storage
Reagent (A): Diff-Quik I	100mL	RT, avoid light
Reagent (B): Diff-Quik II	100mL	RT, avoid light

### Introduction

Diff-Quick Staining is a fast staining method based on Wright Staining, which is one of the common staining methods in cytology. The dye solution is prepared by the rapid staining method recommended by the World Health Organization (WHO). Similar to Wright stain, it is improved by Romanowsky Stain technology principle. The staining result is also very similar to that of Wright's, but the time required for Diff Quik Staining is very short, and the staining can be completed within 90s.

Diff-Quick Stain Kit (No Fixative) is mainly used for blood cell smear, bone marrow smear, vaginal secretion smear and exfoliated cell smear. Diff-Quick Stain Kit (No Fixative) is very suitable for batch dyeing with clear background and no sediment.

### Protocols(*for reference only*)

1. Prepare blood smear or bone marrow smear by routine method, and fix in methanol for 20s after natural drying.
2. Staining with Diff-Quik I for 5-10s (lifting the section up and down for 2-3 times to evenly distribute the dye solution), take it out immediately.
3. Staining with Diff-Quik II for 10-20s (lifting the section up and down for 2-3 times to evenly distribute the dye solution), take it out immediately.
4. After washing with water, observe under the microscope immediately when it is wet, or observe after dehydration in 95% ethanol for 5s, 100% ethanol (I) for 5s, 100% ethanol (II) for 30s, transparency by xylene for 1min and sealing.

### Result

Nucleus, Leukocyte	Deep Blue
Matrix, Lymphocyte	Purple
Cytoplasm, Red Blood Cell	Pink Red

### Note

1. Blood smear or bone marrow smear should be uniform in thickness, so as not to affect the staining effect.
2. Blood cell smear staining requires fresh whole blood or EDTA anticoagulant.
3. Bone marrow smear staining requires the prepared smear to be shaken or dried quickly in the air, so as to prevent cell shrinkage and deformation or hemolysis due to air humidity. It is not allowed to dry the smear with high temperature or fire baking.
4. Smear staining of vaginal secretions requires fresh specimens to be fixed with flame or alcohol as soon as possible after smearing to avoid cell deformation.
5. After natural drying, fix it with dye solution to prevent flaking.
6. Please do not remove the dye solution or wash the smear directly. Do not pour out the dye first, so as to prevent the dye from settling on the smear.
7. The dye solution can be reused, but it can not be repeated many times. If there is sediment, should filter it before use.
8. If the dye is too deep, it can be decolorized properly with methanol or alcohol, and it is better not to be re dyed.
9. If the dyeing is too deep or too shallow, adjust the dyeing time or working solution concentration.
10. The pH value has certain influence on dyeing. The section should be clean and free of acid and alkali pollution, so as not to affect the dyeing effect.