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Rakoff Stain Solution, For Exfoliated Cells

Cat: G1555

Size: 100mL/500mL

Storage: RT, avoid light, valid for 6 months.

Introduction

Conventional cytological stain method includes Papanicolaou Stain, Rakoff Stain and Giemsa Stain. Papanicolaou Stain and Rakoff Stain can be used to detect the estrogen level of vaginal epithelium. Rakoff stain is easy to operate and is often used for rapid determination of estrogen level in clinical vaginal cytology smears.

Protocol(for reference only)

1. Collect the squamous epithelial cells of the mucosa along the side wall of the vagina with cell brush into a test tube containing 1-2 mL physiological saline.
2. Add 3-5 drops of Rakoff Stain Solution into the test tube and mix gently.
3. Add 1-2 drops of mixed solution on the section and make a cell smear.
4. When the smear is to be dried , seal it with resinene.

Result

Acidophilic cytoplasm	Pink
Basophilic cytoplasm	Light green
Keratinocyte pyknosis	No color

Note

1. When making cell smears, pay attention to make cells well-distributed.
2. For your safety and health, please wear experimental clothes and disposable gloves.