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**Cat:** G1650

**Size:** 100mL

**Storage:** RT, avoid light, valid for 6 months.

### Introduction

Methyl green, also known as double green SF, is an alkaline dye. It is a green micro crystal or powder with metallic luster. Its molecular weight is 608.78, and its formula is  $C_{27}H_{35}Cl_4N_3Zn$ . Methyl Green Stain Solution, 0.1% is used to stain the nucleus in tissue or cell staining. It is often used to identify DNA. The DNA in the nucleus will be dyed blue-green when exposed to methyl green. It can also be used for immunofluorescence or immunohistochemical staining. 10 ml staining solution can stain 20 samples.

### Self Provided Materials

4% PFA, Distilled water, Series of ethanol.

### Protocol (for reference only)

#### Sample treatment

For paraffin section, dewax in xylene twice for 5-10mins each. Rehydrate with series of ethanol (100%,90%,80%,70%) for each level 3mins,and finally in distilled water 3min.

For frozen section, restore to RT and wash with distilled water for 2mins.

For cultured cell,fix with 4% PFA for more than 10mins.Then wash with distilled water twice for 2mins each.

#### Methyl green staining

1. For the above treated samples, stain with Methyl Green Stain Solution,0.1% for 5-10mins.
2. Rinse with distilled water twice for 1 min each.Then watch under microscope or seal following the step.
3. (optional) if immunofluorescence staining is required, buffer solution (PBS or saline or TBS, etc.) of diluted fluorescent dye is recommended for 5min.
4. (optional) Immunofluorescence staining was performed as required.
5. Series of ethanol dehydration, 1 min per gradient, xylene transparent twice for 1 min each. Neutral gum seal.

### Result

Nucleus	Green or Blue Green
---------	---------------------

### Note

1. When using the dye solution for the first time, it is recommended to take 1-2 samples for pre test before use.
2. Methyl green staining can adjust the time according to the staining results and requirements.
3. For your safety and health, please wear experimental clothes and disposable gloves.