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Cat: G1995

Size: 3×20mL

Storage: RT, avoid light, valid for 6 months.

Kit components

Reagent	Size	Storage
Reagent (A): Auramine-Phenol Staining Solution	20mL	RT, avoid light
Reagent (B): Differentiation Solution	20mL	RT
Reagent (C): Oxidant Solution	20mL	RT, avoid light

Introduction

Cryptosporidium oocysts staining methods are Auramine-Phenol Staining, Modified Acid Fast Staining and so on. Auramine-Phenol Stain Kit is composed of auramine solution, Differentiation Solution and Oxidant Solution. It is suitable for fresh fecal samples or samples fixed by formalin.

The reagent is only used in scientific research and is not suitable for clinical diagnosis or other purposes.

Protocols(for reference only)

1. Prepare thin smears of fecal specimens, dry them naturally, and fix them in methanol for 2-3min.
2. Drop the Auramine-Phenol Staining Solution on the dried fecal film, place it for 10-15min, and wash it with distilled water.
3. Add The Differentiation Solution onto the fecal membrane, place for 1 min and wash with distilled water.
4. Drop Oxidant on the fecal membrane, place it for 1min and wash it with distilled water.
5. Dry and view under the fluorescence microscope.

Result

Under the low power fluorescence microscope, the egg capsule is milky white or slightly green. The egg capsule wall is a thin layer. Most of the egg capsules are deeply stained around and lightly stained in the center.

Note

1. Nonspecific fluorescent particles can appear in some specimens, which should be distinguished.
2. For your safety and health, please wear experimental clothes and disposable gloves.