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Cat: G2300

Size: 500mL

Storage: RT, avoid light, valid for 6 months.

### Introduction

The fixative can be divided into aldehyde fixative, mercury fixative, alcohol fixative, oxidant fixative, picric acid salt fixative, etc. Formaldehyde in aldehydes and alcohol in alcohols are more commonly used.

BS Fixative is mainly composed of formaldehyde, zinc sulfate and deionized water, which is suitable for tissue fixation of most immunohistochemical staining.

### Protocol(for reference only)

1. Get fresh samples, use saline or 1×PBS was washed and immersed in BS Fixative quickly.
2. The fixation time of the specimens is controlled at 1 ~ 4h/mm, and the fixation time of large specimens should be extended appropriately.
3. Carry out subsequent dewatering and embedding operation.

### Note

1. The fixing time changes due to different thickness of tissue samples. Proper material selection for the tissue is beneficial to the permeation of the fixative.
2. The capacity of the fixative should be sufficient. Generally, the volume ratio of fixative to tissue block should be more than 10:1.
3. After fresh tissue is taken out, it should be fixed in time. If it can not be fixed in time, it should be stored in physiological saline for inspection in time.
4. For your safety and health, please wear experimental clothes and disposable gloves.