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**Cat:** G2542

**Size:** 2×50mL

**Storage:** 2-8°C, avoid light, valid for 1 year.

### Kit Components

Reagent	2×50mL	Storage
Reagent (A): Alcian Acidic Solution	50mL	RT
Reagent (B): Alcian Staining Solution(pH2.5)	50mL	2-8°C, avoid light

### Introduction

Alcian also known as Alcian Blue, is a kind of copper titanium cyan conjugate dye, which was originally used for textile fiber dyeing. This kind of cationic dye combines with acid group, that is, alcian blue forms insoluble complex with anionic group such as carboxyl group and sulfate group. Alcian Blue is composed of the central copper containing phthalocyanine ring and four isothioureia groups linked by thioether bond. The isothioureia group is moderately alkaline, making Alcian Blue with cation. When the pH value is 2.5, the carboxyl group in the tissue is ionized with a negative charge and forms a salt bond with the cation in Alcian Blue, which makes the acid mucus (sulfate mucin and salivary mucin) with carboxyl group stained. When the pH value is 1.0, the sulfate radical in the tissue is ionized with a negative charge, forming a salt bond with the cation in Alcian Blue, so that the tissue with sulfate radical (such as sulfuric acid slime) is stained.

Alcian Blue Cartilage Stain Kit, pH 2.5 can distinguish the genus of mucus by different pH value of staining solution. When pH value is 1, carboxyl group(COOH) is not dyed, and sulfuric acid group(OSO<sub>3</sub>H) is dyed. When pH value is 2.5, COOH is dyed well, while OSO<sub>3</sub>H is dyed bad, which makes mucin sulfate and mucin saliva stain. This reagent is only suitable for scientific research, not for clinical diagnosis or other purposes.

### Self Provided Materials

Series of ethanol, distilled water

### Protocols(for reference only)

1. Conventional decalcification, dewax to distilled water.
2. Soak in Alcian Acidic Solution for 3min.
3. Stain with Alcian Staining Solution(pH2.5) for 30min.
4. Rinse in running water for 5min.
5. Dehydrate in series of ethanol, transparent in xylene and seal with resinene.

### Result

Sulfated Mucin	Blue
Nucleus	Light Blue

### Note

1. The fixative adopts 10% neutral formalin.
2. To selectively identify sulfated mucins and proteoglycans, Alcian Blue with low pH value (pH=1.0) should be used, it is to adjust the pH value to 1.0. The dyeing procedure is the same as that of Alcian Blue with pH =2.5, and the dyeing time should be prolonged accordingly.
3. The opened reagent shall be used up within 6 months after open, and the cap shall be tightened timely after each use to avoid volatilization or deterioration.
4. For your safety and health, please wear experimental clothes and disposable gloves.