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# Endogenous Avidin-Biotin Blocking Kit

**Cat:** G2790

**Size:** 2×100mL

**Storage:** 2-8°C, avoid light, valid for 6 months.

## Kit components

| Reagent  | Size  | Storage |
|--|-------|---------|
| Reagent(A):Endogenous Avidin Blocking Solution | 100mL | 2-8°C   |
| Reagent(B):Endogenous Biotin Blocking Solution | 100mL | 2-8°C   |

## Introduction

Endogenous biotin is widely distributed in mammalian tissues, especially in liver, lung, spleen, fat, brain and other tissues. After formalin fixation, the endogenous background will be reduced, but the background in frozen sections will be greatly enhanced. Intense antigen repair will expose endogenous substances in formalin fixed tissues, such as endogenous avidin-biotin. The kit contains endogenous avidin blocking solution and biotin blocking solution, which can well inhibit the activity of endogenous substances (mainly endogenous avidin-biotin), especially suitable for blocking the activity of endogenous substances in immunohistochemistry.

## Self Provided Materials

PBS, Dyeing Rack, Incubator.

## Protocols(*for reference only*)

1. Slice dewaxing and hydration.
2. Add Endogenous Biotin Blocking Solution and incubate at 37 °C for 15-30min.
3. Fully wash with PBS buffer for 5min.
4. Add Endogenous Avidin Blocking Solution and incubate at 37 °C for 15-30min.
5. Fully wash with PBS buffer for 5min.
6. Conduct steps of immunostaining.

## Note

1. If the tissue to be tested is rich in endogenous streptavidin biotin, the non streptavidin biotin detection system should be considered.
2. Do not use the reagent if there is more turbidity.
3. Avoid repeated freezing and thawing, otherwise it is easy to fail.