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规格: 100mL/500mL

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**Cat:** G3004

**Size:** 100mL/500mL

**Storage:** 2-8°C, avoid light, valid for 1 year.

### Introduction

2,3,5-triphenyltetrazolium chloride (TTC) is a lipid soluble light sensitive complex, which is commonly used to detect the viability of seeds and the ischemic infarction of mammalian tissues. It is the proton receptor of pyridyl-nucleoside structural enzyme system in respiratory chain, which can react with dehydrogenase in normal tissues and show red color. The dehydrogenase in the living cells of the viable seed embryo tissue can reduce TTC to insoluble red stable triphenylmethane (TTF). If the embryo dies or the viability of the embryo declines, it cannot be stained or stained shallowly. Therefore, the viability of seeds can be determined according to the location and depth of staining.

### Protocols(*for reference only*)

1. Soak the seeds to be tested in warm water (30-35 °C) for 2-6h to make the seeds fully imbibed.
2. Take 100 seeds randomly, cut them into two halves along the center line of embryo, and put half of each seed in two petri dishes for use.
3. (*optional*)place half of the seeds in boiling water and boil for 5min to kill cells as negative control.
4. Add TTC Solution, 0.4% to the petri dish containing seeds to submerge the seeds.
5. Put it into a constant temperature box at 30-35 °C for 30min-1h.
6. Pour out TTC Solution, 0.4%, wash with tap water for 2-3 times, and immediately observe the staining of seeds.

### Result

Active seeds	Red
Negative control and inactive seeds	Not colored

### Note

1. TTC Solution, 0.4% can not be reused if it turns red.
2. The seeds should be absolutely immersed in TTC Solution, 0.4% when staining.
3. Immediately observe after staining with TTC Solution, 0.4% and the seeds can not be placed for a long time.
4. For your safety and health, please wear experimental clothes and disposable gloves.