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## Luxol Fast Blue Myelin Stain Kit(Cresyl Violet Method)

**Cat:** G3245

**Size:** 3×50mL

**Storage:** RT, avoid light, valid for 1year.

### Kit Components

Reagent	3×50mL	Storage
Reagent(A): LFB Staining Solution	50mL	RT, avoid light
Reagent(B): LFB Differentiation Solution	50mL	RT
Reagent(C): Cresyl Violet Staining Solution	50mL	RT, avoid light

### Introduction

Myelin sheath is a multilayer lipid double-layer structure formed by the plasma membrane of nerve cells spirally winding along the axis of axon. Myelin staining has certain significance in pathological diagnosis. The pathological changes of myelin sheath can be divided into early, middle and late stages. In the early stage of the lesion, the color is deep. In the middle stage, the myelin degenerated into lipid droplets. In the late stage, the myelin sheath degenerated completely and is cleared by phagocytes, so there is no positive result.

Many diseases can cause the change of myelin sheath. Luxol Fast Blue Myelin Stain Kit can show whether the myelin sheath is complete, denatured, necrotic and repaired under pathological conditions. It has significance for the pathological diagnosis and research of nerve tissue.

### Protocols(for reference only)

1. Cut into paraffin section in 5-8 $\mu$ m thick, then dewax to 95% ethanol.
2. Soak in LFB Staining Solution and stay overnight at room temperature (the staining time of frozen section shall not exceed 16h). Wash the excess dye solution with 95% ethanol. Rinse in distilled water.
3. Differentiate by LFB Differentiation Solution for 15s. Differentiate in 70% ethanol for 30s.
4. Rinse in distilled water and observe under the microscope(if necessary, can repeat step 3)until the gray and white matter is clear.
5. Re-dyeing with Cresyl Violet Staining Solution for 30-40s, wash with water.
6. Dehydrate in 95% and 100% ethanol, transparent by xylene and seal with resinene.

### Result

Myelin Sheath	Blue
Neuron	Pink Red to Purple

### Note

1. Differentiation is a key step. The differentiation time should be strictly controlled, and the degree of differentiation can be observed under the microscope.
2. It is better to use 10% formalin as fixative.
3. The section should be controlled within 8-9  $\mu$ m, otherwise, it is easy to take off the section or over dye.
4. If the dyeing effect at room temperature is not good, can stain at 56 °C for 2h.
5. After washing away the re-dyeing solution, it is not advised to dehydrate in 70% ethanol, otherwise the blue color of Luxol fast blue will be removed.
6. For your safety and health, please wear experimental clothes and disposable gloves.