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Myelin Stain Kit (Fast Green FCF Method)

Cat: G3246

Size: 3×50mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent	3×50mL	Storage
Reagent (A): Fixative	100mL	RT, avoid light
Reagent (B): Fast Green FCF Staining Solution	50mL	RT, avoid light
Reagent (C): MS Differentiation Solution	100mL	RT

Introduction

Myelin Sheath is a layer of membrane wrapped around the axon of nerve cells, that is, myelin sheath is composed of myelin cells and cell membrane. It is a multilayer lipid double-layer structure formed by the plasma membrane of nerve membrane cells spirally winding along the axis of axon. There is a Langfei's node on myelin sheath, which can make nerve impulse jump and transmit. Myelin staining has certain significance in pathological diagnosis. The pathological changes of myelin can be divided into early, middle and late stages. In the early stage, the color is deep. In the middle stage, the myelin degenerates into lipid droplets, which can be displayed by lipid staining. In the later stage, the myelin degenerates completely and is removed by phagocytes, so there is no positive result of myelin sheath.

Many diseases can cause the change of myelin sheath. Myelin Stain Kit (Fast Green FCF Method) can show whether the myelin sheath is complete, degenerated and also can show necrotic degree and repair situation under pathological conditions. It has significance for the pathological diagnosis and research of nerve tissue. The myelin sheath is dark green and the demyelinated fiber is not stained.

Protocol (for reference only)

Fixation: fix the samples with formaldehyde calcium fixative.

1. For paraffin section, dewax to distilled water.
2. Slightly wash with 95% ethanol for 1 min.
3. Put the section in Fast Green FCF Staining Solution and dye in 37 °C incubator for 30-45min.
4. Directly wash with 95% ethanol twice for 10s each time.
5. Rinse with distilled water.
6. Differentiate by MS Differentiation Solution for 10s.
7. Rinse with distilled water (repeat step 6 if differentiation is insufficient).
8. Conventionally dehydrate in series of ethanol and transparent by xylene, seal with resinene.

Result

Myelin Sheath	Dark green, annular in cross section, cord or fishbone in longitudinal section.
Demyelinated Fiber	No color, semi-circular or non color in cross section, show blank area.

Note

1. Differentiation is a key step. The differentiation time should be strictly controlled, and the degree of differentiation can be observed under the microscope.
2. It is better to use the fixative provided by this kit for the specimen. If it is self prepared, 10% formaldehyde calcium fixative can be used.
3. The section should not be too thick, and should be controlled within 5-7 μ m, otherwise, it is easy to take off the section or over dyeing.
4. For your safety and health, please wear experimental clothes and disposable gloves.