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## Gomori Ammoniacal Silver Stain Solution

**Cat:** G3530

**Size:** 100mL

**Storage:** 2-8°C, avoid light, valid for 6 months.

### Introduction

Reticular fiber is a kind of fiber in reticular connective tissue, which is produced by reticular cells. Its diameter is mostly 0.2-1.0 $\mu$ m. It has toughness but no elasticity. There are many dyeing methods for reticular fiber, but the dyeing principle is basically the same, and most of them adopt ammonia silver staining method. The dyeing principle of Gordon-Sweets staining is that ammonia silver solution is easy to be adsorbed by tissue and combined with protein of tissue. It is reduced to black or brown black metallic silver by formaldehyde and deposited in tissue and on its surface. In the traditional method, after reduction, gold chloride solution is used for color matching, and then hypo solution is used to wash away the unreduced silver salt on the tissue. This improved method omits this step and makes the contrast of reticular fiber clearer.

Gomori Ammoniacal Silver Stain Solution is the main dyeing component of G3535-Reticulin Stain Kit(Gomori Silver Method). Frozen section, low temperature section and collodion section can be used for reticular fiber staining. All kinds of fixatives can be used, and the heavy metal mercury salt or osmium salt fixatives occasionally produce some non-specific silver background.

It is often used to distinguish the nature and source of tumor, cancer and sarcoma, lymphosarcoma and reticulosarcoma, hemangioendothelioma and hemangiopericytoma, etc.

### Protocols(for reference only)

1. Fix the tissue in 10% formalin fixative, then conventionally dehydrate and embed.
2. Cut into paraffin section in 4 $\mu$ m thick, and dewax to distilled water.
3. Put the section on the dyeing frame, drop in Gomori Oxidant and oxidize for 5min, wash with water slightly.
4. Bleach in Oxalic Acid Solution for 1-2min, rinse in running water for 2min and slightly with distilled water.
5. Mordant dyeing with Ferric Ammonium Sulfate Solution for 5min, wash slightly with tap water and distilled water.
6. Add Gomori Ammoniacal Silver Stain Solution and dye for 3-5min, wash slightly with distilled water.
7. Reduce by Gomori Reductant for 5-10min, rinse in running water for 10min.
8. Conventionally dehydrate and transparent, seal with resinene.

### Result

Reticular Fiber	Black
Collagen Fiber	Yellow or Yellowish Brown
Nucleus	Brown or Black Brown

### Note

1. Glass container must be soaked in washing solution for one day, washed with tap water and washed twice with distilled water.
2. 10% formalin fixative is a suitable fixative. It is not suitable to use fixative containing mercury such as Zenker's Solution, otherwise it is easy to cause nonspecific deposition of sections.
3. Gomori Ammoniacal Silver Stain Solution is not very stable and sensitive to light. It should be kept away from light at 4 °C and used after restoring to room temperature.
4. For your safety and health, please wear experimental clothes and disposable gloves.