Griess 试剂

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[1] Yu Bao,Zhi Li,Si-Hao Chen,et al.Artemisinin is highly soluble in polyethylene Glycol 4000 and such solution has multiple biological effects.Acta Biochimica Polonica, Quarterly by Polish Biochemical Society.May 2020.(IF 1.626)

Griess Reagent

Cat: G9440

Size: 10g

Storage: RT, avoid light, valid for 6 months.

Introduction

Griess experiment is mostly used for qualitative experiment of urine nitrite. Its principle is that gram negative bacteria such as Escherichia coli can reduce nitrate in urine as nitrite. Nitrite diazotizes p-aminobenzenesulfonic acid in Griess reagent to p-diazobenzenesulfonic acid, and diazopine sulfonic acid is combined with α -naphthylamine to form n- α -naphthylaminobenzenesulfonic acid, which is red. The reagent is only used in scientific research and is not suitable for clinical diagnosis or other purposes.

Protocol(for reference only)

1. Add 5mL of urine into the test tube and add 0.05g Griess Reagent immediately.
2. Oscillation mixing, observation results.

Result

Orange Red Or Dry Powder Turned Red	Beyond the Measurement Limit
Pink To Rose	Positive
No Change In Color	Negative

Note

1. Urine must be fresh.
2. The positive results are determined by the following four conditions:
 - ① There is a proper amount of nitrate in the body;
 - ② The bacteria in urine contain nitrate reductase;
 - ③ The retention time of urine in bladder should be more than 4 h;
 - ④ The use of antibiotics may interfere with the experiment within 48 h.

Reference

[1] Yu Bao,Zhi Li,Si-Hao Chen,et al.Artemisinin is highly soluble in polyethylene Glycol 4000 and such solution has multiple biological effects.Acta Biochimica Polonica, Quarterly by Polish Biochemical Society.May 2020.(IF 1.626)