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Lipid Fluorescent Staining Kit (Nile Red Method)

Cat: G1264

Size: 2×50mL

Storage: -20°C, avoid light, valid for 1 year.

Kit Components

| Reagent | | 2×50mL | Storage |
|--|--------------------------|--------|--------------------|
| Reagent(A): Lipid Fixative | | 50mL | 2-8°C, avoid light |
| Reagent(B): Stain Working Solution | B1:Dyeing Stock Solution | 0.2mL | -20°C, avoid light |
| | B2:Buffer Diluent | 40mL | 2-8°C |
| Before the use, prepare the dyeing working liquid according to the ratio of B1: B2=1:200, pay attention to the light avoidance operation, and finish using it within 1h. | | | |

Introduction

Lipids are a natural mixture of low polar (hydrophobic) natural macromolecules, including fatty acids and their derived glycerides and phospholipids, as well as steroids, such as cholesterol. The formation of cytoplasmic lipid droplets is a normal cell physiological process, consisting of neutral lipids, usually triglyceride lipid droplets, as energy reserves; a few are cholesterol ester lipid droplets, as the storage of excessive cellular cholesterol. Nile red (Nile Red) stain is a lipophilic oxazine yellow light dye. In combination with lipids such as wax ester, triglycerides and various fatty acids, a strong orange fluorescence (emission wavelength of 598nm) was displayed at an excitation wavelength of 543nm. Also shows red under irradiation with UV light.

The advantage of the Lipid Fluorescent Staining Kit (Nile Red Method) is to provide a composite fluorescent staining solution for simultaneous tissue in situ staining of lipids and nuclei. Under green light excitation, the orange-colored lipid droplets, and the blue cell nuclei upon UV excitation, were suitable for frozen sections and cell samples. This product is ready to use kit, stable performance, clear and sensitive coloring.

Self Provided Materials

1× PBS, 1.5 ml Centrifuge tube, Fluorescence microscope

Protocols(for reference only)

Before use, take out the reagent (B1) and reagent (B2) and put them back to room temperature. Prepare the dyeing working solution according to the ratio of B1: B2 = 1:200, avoid the light, and then carry out the following operations.

For Frozen Section Stain

1. Take out the prepared frozen section, place it for 3min, and rewarm it to room temperature.
2. Add the Reagent(A): Lipid Fixative to the section, cover the whole tissue and fix at room temperature for 15-20min. Wash with 1×PBS 3 times, 1min for each.
3. Add the Stain Working Solution to the sections and incubate for 10min at RT avoid light.
4. Wash the section 1-2 times with 1×PBS for 30s each til no red pigment comes out.
5. Cover slides were added to prepare temporary sections or use S2100-antifluorescence decaying tablets, and then observed under a fluorescence microscope: excitation wavelength 543nm, emission wavelength 598nm (lipid droplet staining) and ultraviolet light (nuclear staining).

For Adherent Cells Stain

1. Remove the medium from the dish or six-well plate carefully.
2. Add 1ml 1×PBS carefully, clean the cell surface gently and repeated 1-2 times, each for 30s.

3. Remove the cleaning solution from the culture dish or from the six-well plate carefully.
4. Add 1ml Reagent (A): Fixation Solution, cover the cell surface and fix for 10-15min. Remove the fixative fluid carefully.
5. Add 1ml 1×PBS and clean the cell surface, repeated 1-2 times. Remove the washing liquid carefully.
6. Add appropriate amount of staining solution and incubated for 10min at room temperature in the dark.
7. Add 1ml 1×PBS carefully, and wash the cell surface for 30s-1min by shaking, repeat 1-2 times.
8. Cleaning liquid in the last step is not abandoned, pay attention to avoid light and observe by fluorescence microscope within 30 min: excitation wavelength 543nm, emission wavelength 598nm (lipid droplet coloring) and ultraviolet light (nuclear staining).

Result

| | |
|------------------|-----------------------------------|
| Lipid substances | Red or orange-yellow fluorescence |
| Nucleus | Bright blue fluorescence |

Note

1. The Stain Working Solution should be configured before the use , not in advance.
2. Depending on the tissue and cell differences, the dye dilution ratio (1:100-1:400) can be adjusted appropriately to enhance or decrease the fluorescence intensity.
3. Avoid using lipids in the cell culture medium. During staining, incubation and observation, avoid intense light irradiation and avoid fluorescence quenching.
4. When you use this kit, the reagent (B1): Dyeing Stock Solution should be avoid repeated freezing and thawing, can be divided into small specifications for preservation .
5. It is recommended that fluorescence detection and analysis be performed immediately after the staining is completed. If there is no fluorescence microscope, UV observation can be selected.
6. For your safety and health, please wear lab suits and wear disposable gloves.