DiI-Low Density Lipoprotein DiI 标记低密度脂蛋白

Low Density Lipoprotein labeled with 1,1'-dioctadecyl-3,3,3',3'-tetramethyl-indocarbocyanine perchlorate

Cat No: IL2120

Concentration: >1mg/mL protein (Based on actual label concentration) .Buffer containing phosphate-buffered saline at pH 7.4 and 0.2 mM EDTA.

Specifications: 0.5mg/vial; 0.22 micron membrane filtered, aseptically filled. Cell Culture Tested.

Absorbance Ratio: DiI/Protein=555nm/275nm=1.40

Storage

This product is stable for 6 weeks after receipt when handled aseptically and stored at 2-8°C (Do not freeze). After prolonged storage, some precipitate may be observed. This is normal for this product. Clarify out the aggregates by spinning in centrifuge tube for 2 minutes. LDL products have a natural tendency to aggregate. Aggregates of this product can interfere with its use.

Product Preparation:

Purified LDL is labeled with the fluorescent probe, DiI, and reisolated by ultracentrifugation (1.019-1.063). Each lot is evaluated on a murine macrophage cell line for fluorescence uptake.

Typical Lipoprotein Labeling Protocol

1. Aseptically dilute the DiI-LDL to 10-40 µg/mL in growth media.
 2. Add to live cells and incubate for 4-6 hours at 37°C.
 3. Remove media.
 4. Wash cells several times with probe-free media.
- A. Fluorescence Microscopy:
Visualize using standard rhodamine excitation: emission filters (or suggested wavelengths excitation:emission at 549nm:565nm). If fixation is desired use 3% formaldehyde in PBS. (Do not use methanol or acetone fixation - DiI is soluble in organic solvents).
Note: A positive culture must be stained for comparison purposes.
- B. Cell Sorting:
Label as in steps 1-5. Trypsinize or treat cultures with EDTA to produce a single cell suspension. Use labeled pure cultures of positive and negative cell types to set gates on the cell sorter.
Suggested wavelengths for cell sorting: Excitation: 514/549nm. Emission: 565nm.

Fixation and Mounting of DiI Labeled Cells

1. Wash 3 times in PBS.
2. Fix in 3% formaldehyde/PBS for 20 minutes at room temperature.
3. Rinse 5 seconds in distilled water at room temperature.
4. Drain liquid onto chem-wipe.
5. Invert cover, slip on a drop of 90% Glycerol and 10% PBS onto a microscope slide.
6. Seal with Kroenigs wax, also known as cover glass cement Do not use nail polish. Store at -20°C

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