DiI-Ox-Low Density Lipoprotein DiI标记氧化型低密度脂蛋白

Cat : IL2100

Size: 500µg protein/vial

Concentration: > 1 mg/mL protein (Based on actual label concentration)

Storage:

This product is stable for 6 weeks after receipt when handled aseptically and stored at 2-8°C (Do not freeze). After prolonged storage, some precipitate may be observed. This is normal for this product. Clarify out the aggregates by spinning in centrifuge tube for 2 minutes.

Introduction:

DiI-Ox-LDL, oxidized Low Density Lipoprotein, labeled with 1, 1'-dioctadecyl-3,3,3',3'-tetramethyl-indocarbocyanine perchlorate, labels both vascular endothelial cells and macrophages. It can be used to identify and/or isolate these cells from mixed cell populations and investigate uptake of modified LDL by different cell types. When cells are labeled with DiI-Ox-LDL, the lipoprotein is degraded by lysosomal enzymes and the DiI (fluorescent probe) accumulates in the intracellular membranes. Labeling cells with DiI-Ox-LDL has no effect on cell viability. Pure cultures of vascular endothelial cells can be isolated from complex primary cultures using fluorescent activated cell sorting based on their increased metabolism of the DiI-Ox-LDL. Contaminating cell types (fibroblasts, smooth muscle, pericytes, epithelial cells) are not labeled. Macrophages can be differentiated from mixed cell populations (including endothelial cells) because they are more brightly labeled. Labeling endothelial cells with DiI-Ox-LDL has many advantages over labeling other endothelial cell associated antigens. The labeling procedure is one step, and once the cells are labeled, the fluorescent probe (DiI) is not removed by Trypsin. Both low density and confluent cultures of vascular endothelial cells are effectively labeled. No other cell type (other than macrophages) is labeled to the same level as vascular endothelial cells. Each lot of DiI-Ox-LDL is evaluated for the specific labeling of bovine aortic endothelial cells and murine macrophages to assure consistent results.

Protocol

1. Aseptically dilute the DiI-OxLDL to 20-40 µg/mL in your culture media.
2. Add to live cells and incubate for 3-6 hours at 37°C.
3. Remove media containing DiI-OxLDL from your culture.
4. Wash cells several times with probe-free media.
5. A. Fluorescence Microscopy:

Visualize using standard rhodamine excitation: emission filters (or suggested wavelengths excitation:emission at 549nm:565nm). If fixation is desired use 3% formaldehyde in PBS. (Never use methanol or acetone fixation - DiI is soluble in organic solvents).

Note:

A positive culture must be stained for comparison purposes.

B. Cell Sorting:

Label as in steps 1-5. Trypsinize or treat cultures with EDTA to produce a single cell suspension.

Use labeled pure cultures of positive and negative cell types to set gates on the cell sorter.

Suggested Wavelengths for Cell Sorting: Excitation: 514/549nm; Emission: 565nm

Fixation and Mounting of DiI Labeled Cells

1. Wash 3 times in PBS.
2. Fix in 3% formaldehyde/PBS for 20 minutes at room temperature.
3. Rinse 5 seconds in distilled water at room temperature.
4. Drain liquid onto chem-wipe.
5. Invert cover, slip on a drop of 90% Glycerol and 10% PBS onto a microscope slide.
6. Seal with Kroenigs wax. Do not use nail polish. Store at -20°C.

Special Note:

- * LDL products have a natural tendency to aggregate. Aggregates of this product can interfere with its use. To clarify these aggregates out, simply spin in a microfuge for 2 minutes.
- * Product information is for reference only. Welcome to call 400-968-6088 for consultation.
- * This product is for scientific research use only. Not used for medical, clinical and other purposes.
- * For your safety and health, please wear lab coats and wear disposable gloves and masks

References:

(See more information on <http://www.solarbio.com/>)

(Welcome to follow our company scholarship program: <http://www.solarbio.com/lw.php>)

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