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**Storage:** 2-8°C, valid for 1 year.

### Introduction

This product is a kind of dye used to prepare SDS-PAGE stacking gel. The dye is available in four colors: red, green, black and blue.

The addition of this dye can make the sample holes of gel protein concentrate clear at a glance, facilitate the accurate addition of protein samples into the sample holes, and make it easy to judge the edge distance of the loading holes or whether they are distorted or damaged. And the color of the dye will not enter the separating gel with electrophoresis, and it can stably exist in the stacking gel. After electrophoresis, it is also easy to identify and cut off the stacking gel, without affecting the electrophoresis and dyeing effects. It can be used with our P1200-SDS-PAGE Gel Kit.

### Protocols(for reference only)

*Due to the special nature of the dye, precipitation will occur during storage, so please turn upside down or blow it to mix well before use, and then it can be used normally.*

1. Taking the 1.0 mm thick gel as an example, it takes about 2 ml of stacking gel solution to prepare a single piece of gel. As the ratio of 1:200, add 10  $\mu$ L Dye for PAGE Stacking Gel, 200 $\times$  into 2ml stacking gel solution and mix evenly (the amount of pigment added can be determined according to personal preferences, and the maximum amount of pigment used should not exceed 5 times of the recommended amount).
2. Add gel polymerization agent (such as ammonium persulfate or other ammonium persulfate substitutes) and catalyst (TEMED or TEMED substitutes). After fully mixing, use a pipette to inject the stacking gel solution between the glue making glass plates and insert an electrophoresis comb.
3. After the stacking gel solidifies, pull out the electrophoresis comb for subsequent electrophoresis.

### Note

1. The dye can also be added after adding the gel polymerization agent (such as ammonium persulfate or other ammonium persulfate substitutes) and catalyst (TEMED or TEMED substitutes), and then directly pour the stacking gel after mixing.
2. Since it is common for different samples to electrophoresis multiple gel at the same time, we provide stacking gel dyes of multiple colors, including red, green, black and blue, to distinguish different samples or different gel electrophoresis of different operators.
3. This product is only for scientific research of professionals, and shall not be used for clinical diagnosis or treatment, food or drugs, or stored in ordinary houses.
4. For your safety and health, please wear experimental clothes and disposable gloves