DiO-Low Density Lipoprotein DiO 标记低密度脂蛋白

Cat No: IL2190

Concentration: >1mg/mL protein (Based on actual label concentration) .Buffer containing phosphate-buffered saline at pH 7.4 and 0.2 mM EDTA.

Specifications: 0.5mg/vial; 0.22 micron membrane filtered, aseptically filled. Cell Culture Tested.

Excitation:484nm

Emission:506nm

Storage

This product is stable for 6 weeks after receipt when handled aseptically and stored at 2-8°C (Do not freeze). After prolonged storage, some precipitate may be observed. This is normal for this product. Clarify out the aggregates by spinning in centrifuge tube for 2 minutes. LDL products have a natural tendency to aggregate. Aggregates of this product can interfere with its use.

Typical Lipoprotein Labeling Protocol

1. Aseptically dilute the Dio-LDL to 10-40 ug/mL in growth media.
2. Add to live cells and incubate for 4-6 hours at 37°C
3. Remove media.
4. Wash cells several times with probe-free media.

A. Fluorescence Microscopy:

Visualize using standard rhodamine excitation: emission filters (or suggested wavelengths excitation :emission at :484nm:506nm). If fixation is desired use 3% formaldehyde in PBS. (Do not use methanol or acetone fixation - DiO is soluble in organic solvents).

Note: A positive culture must be stained for comparison purposes.

B. Cell Sorting:Label as in steps1-5. Trypsinize or treat cultures with EDTA to produce a single cell suspension. Use labeled pure cultures of positive and negative cell types to set gates on the cell sorter.

Suggested wavelengths for cell sorting: Excitation: 484nm. Emission: 506nm.

Fixation and Mounting of DiO Labeled Cells

1. Wash 3 times in PBS.
2. Fix in 3% formaldehyde/PBS for 20 minutes at room temperature.
3. Rinse 5 seconds in distilled water at room temperature.
4. Drain liquid onto chem-wipe.
5. Invert cover, slip on a drop of 90% Glycerol and 10% PBS onto a microscope slide.
6. Seal with Kroenigs wax, also known as cover glass cement. Do not use nail polish. Store at -20°C

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