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Collagen Fiber And Elastic Fiber Staining Kit(EVG-Weigert Method)

Cat: G1590

Size: 3×50mL/3×100mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit Components

Reagent	3×50mL	3×100mL	Storage
Reagent(A): Weigert Fuchsin Staining Solution	50mL	100mL	2-8°C, avoid light
Reagent(B): Acid Differentiation Solution	50mL	100mL	RT
Reagent(C): Modified VG Staining Solution	50mL	100mL	RT, avoid light

Introduction

Elastic fibers are mainly distributed in the arterial walls, alveolar walls, and skin of the human body. When fresh, they are yellow and have strong refractive properties. In routine tissue staining, it is often necessary to simultaneously display two types of fibers to study the coexistence relationship between fibers and their dependence on other structures. Based on this, an EVG composite staining kit was developed on the basis of a separate fiber staining kit, which is used in combination with conventional elastic fiber staining and VG collagen staining.

This reagent kit adopts a combined staining mechanism of weigert fuchsin and modified VG staining solution, which can simultaneously stain elastic fibers and collagen fibers on a single slice without fading easily. After staining, the elastic fibers are purple black, the collagen fibers are red, and the muscles are yellow. It is suitable for observing whether there are abnormal lesions in the elastic fibers and collagen fibers in the tissue, thus assisting in diagnosis. This reagent kit has excellent quality, simple operation, stable staining, and beautiful and clear colors.

Protocol(for reference only)

1. Take fresh samples, fix and embed with paraffin. Cut the section into 3-8um.
2. After 50min of heat treatment on paraffin slices, dewax by xylene for 20min, then treat with gradient ethanol to distilled water for each step 5min.
3. Add into Weigert Fuchsin Staining Solution, cover and dip at room temperature for 1-3 hours, or stain at 56 °C for 1 hour.(see note 2)
4. Differentiate with Acid Differentiation Solution until no staining solution is removed. Wash with distilled water to stop differentiation.(see note 3)
5. Dye with Modified VG Staining Solution for 5-10min, and wash off excess staining solution with distilled water.(see note 4)
6. 75%, 85%, 95%, and 100% ethanol were subjected to conventional dehydration for 3-5 seconds each. Xylene was transparent twice for 1min each time, and neutral gum was used for sealing.
7. Microscopic observation of staining results.

Result

Elastic fiber	Brown or Purple Black
Collagen fiber	Red
Muscle fiber、Red blood and Background	Yellow

Note

1. Most fixing solutions are suitable, but fixing solutions containing chromium salts can result in lighter staining and easier dispersion.
2. If the Weigert Fuchsin Staining Solution is not heated, the staining time can be appropriately extended.
3. The Acid Differentiation Solution usually takes only 3-10 seconds to remove the background color, and the differentiation must be terminated as soon as possible when the color is clear through visual inspection.
4. The Modified VG Staining Solution is stable in dyeing, not easy to fade after washing, and can be thoroughly cleaned.
5. For your safety and health, please wear laboratory clothes and disposable gloves for operation.