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**Cat:** G2210

**Size:** 500mL

**Storage:** 2-8°C, avoid light, valid for 6 months.

### Kit Components

Reagent	500mL	Storage
Reagent (A): Oxidizer Stock Solution	100mL	2-8°C, avoid light
Reagent (B): Dilution Buffer	400mL	2-8°C
Before use, dilute the buffer solution by mixing it in a ratio of A: B=1:4, and use it immediately.		

### Introduction

The fixative can be divided into aldehyde fixative, mercury fixative, alcohol fixative, oxidant fixative, picric acid salt fixative, etc. Formaldehyde in aldehydes and alcohol in alcohols are more commonly used.

PLPD fixative is a lysine buffered oxidative fixative, mainly composed of lysine buffer and oxidative fixative components. It is mainly used for the fixation of carbohydrate containing tissues, especially for the protection of antigenicity and ultrastructure in immune microscopy.

### Protocol(for reference only)

1. Operate according to the specific requirements of the experiment.
2. Generally fix at 2-8°C for 6-18 h.
3. If the tissue block is large, the fixing time can be extended appropriately, but it should not over 24 h.

### Note

1. PLPD Fixative is harmful for human body, please operate carefully in ventilated environment.
2. The oxidant reserve solution is unstable, and it is not recommended to use it if it becomes cloudy or green.
3. If it is not used for a long time (interval of more than 1 month). It is recommended to store it in a cool and dark place.
4. The fixing time changes due to different thickness of tissue samples.
5. The capacity of the fixative should be sufficient. Generally, the volume ratio of fixative to tissue block should be more than 10:1.
6. The effect of temperature on fixation is obvious. Increasing the temperature can accelerate the fixation, but the temperature should not be too high.