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Oil Red O Stain Kit, For Cultured Cells

Cat: G1262

Size: 4×20mL/4×50mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit Components

Reagent		4×20mL	4×50mL	Storage
Reagent(A): ORO Fixative		20mL	50 mL	RT, avoid light
Reagent(B): ORO Stain	B1: ORO Stain A	12 mL	30 mL	2-8°C, avoid light
	B2: ORO Stain B	8 mL	20 mL	RT
Mix B1 with B2 as the ratio of 3:2 and place for 10 min to form ORO Stain, filter it before use.				
Reagent(C): Mayer Hematoxylin Solution		20mL	50 mL	2-8°C, avoid light
Reagent(D): ORO Buffer		20mL	50 mL	RT

Introduction

Neutral fat stains often use Sudan II, Sudan III, Sudan IV, Sudan black B, Oil red O and so on. Sudan dye is often used in traditional methods. Recently, azo dye oil red O is more suitable for dyeing fat. Oil Red O is a strong lipid solvent and dye. It is easy to bind to triglycerides in droplet shape, but has a slightly poor binding to phospholipids. The dyeing principle is generally regarded as a physical miscibility or adsorption, and the fat is dyed by miscibility. The solubility of dyes in frozen sections is higher than that in the original solvent, so when dyeing, the dyes are transferred from organic solvents to lipids and the fats are dyed.

Oil Red O Stain Kit, For Cultured Cells is shorten as ORO Stain and can show least fatty droplet. The sample should not fixed with ethanol (10% formalin if required). The positive staining result of fat is orange to red, but the specific color depended on the concentration of lipid.

Self Provided Materials

60% isopropanol, Distilled water

Protocol(for reference only)

For Cultured Cells

1. Remove the cell culture medium, wash it twice with PBS, and fix it with ORO Fixativ for 20-30min.
2. Discard the ORO Fixative and wash it twice with distilled water.
3. Add 60% isopropanol and soak for 20-30s.
4. After discarding 60% isopropanol, add the newly prepared oil red O dyeing solution and dip it for 10-20min.
5. Discard the staining solution and rinse with 60% isopropanol for 20-30s until the stroma is clear. Wash with water 2-5 times until there is no excess dye.
6. Add Mayer hematoxylin staining solution and counterstain the nucleus for 1-2min. After discarding the dye solution, wash it with water for 2-5 times. Add oil red O buffer for 1min and discard it.
7. The cells were covered with distilled water and observed under a microscope.

For Cell Smear

1. Prepare fresh bone marrow and blood smear, then fix in ORO Fixative for 10-15 min.
2. Take out the smear and put it in the air for 10-15 min.
3. Soak the smear in ORO Stain and dye for 15 min.
4. Rinse with 60% isopropanol for 20-30s, rinse with running water, and wash slightly with distilled water. Re-dyeing the nucleus with Mayer Hematoxylin Solution for 1-2 min.
5. Soak in ORO Buffer for 1min and observed under a microscope..

Result

Neutral Fat	Orange or Red
Phospholipid	Pink Red
Nucleus	Blue

Note

1. ORO Stain is not stable enough and easy to precipitate, so it is not suitable to prepare it in advance.
2. Mayer's Hematoxylin Solution should not stain for too long.
3. The dyeing results can not last for a long time, and should be viewed and photographed as soon as possible.
4. For your safety and health, please wear experimental clothes and disposable gloves.