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Tartrate-Resistant Acid Phosphatase (TRAP) Stain Kit

Cat: G1492

Size: 4×10mL/4×20mL

Storage: -20°C, avoid light, valid for 3 months.

Kit Components

Reagent	4×10mL	4×20mL	Storage	
Reagent(A): TRAP Fixative	50mL	100mL	2-8°C, avoid light	
Reagent(B): TRAP Incubation Solution	B1: AS-BI Buffer	1mL	2×1mL	-20°C, avoid light
	B2: GBC Solution	0.1mL	0.2mL	-20°C, avoid light
	B3: TRAP Buffer	9mL	18mL	RT, avoid light
Mix reagent B1,B2,B3 in 10:1:90 ratio as TRAP Incubation Solution before use and not store for long.				
Reagent(C): Hematoxylin Solution	10mL	20mL	2-8°C, avoid light	
Reagent(D): Methyl Green Solution	10mL	20mL	RT, avoid light	

Introduction

Acid phosphatase (ACP) is widely distributed in various tissues, mainly in the lysosome, so it is often used as a lysosome marker enzyme. The acid phosphatase outside of lysosome exists in endoplasmic reticulum and cytoplasm. Acidic phosphatase is different in all kinds of animals. The optimum pH for its activity is 4.5-5.5. Tartrate-resistant acid phosphatase (TRAP), which mainly exists in normal human alveolar macrophages and leukemia human spleen cells, is released into the cell follicles rather than into the blood. The majority of TRAP in blood comes from osteoclasts, so we can know the function of osteoclasts by measuring TRAP in blood.

The reaction principle of Tartrate-Resistant Acid Phosphatase (TRAP) Stain Kit is that phosphoric acid and naphthol are released by hydrolysis of acid phosphatase with AS-BI as substrate at acidic pH. Naphthol was coupled with diazo salts to form colored products, which were localized in the cytoplasm. If intracellular ACP has tartaric acid-resistant activity, it is positive. It usually used for fresh blood smear, cell smear, frozen section, etc.

Self Provided Materials

Distilled Water, Incubator, Slide, Microscope

Protocol (for reference only)

For Blood or Cell Smear

1. Dry blood or cell smear in air then fix in TRAP Fixative at 4 °C for 30s-3mins and usually for 30-60s.
2. Wash with distilled water and slightly dry in air.
3. Add the TRAP Incubation Solution to sections and place in 37°C incubator and incubate in dark for 45-60min. Then wash with distilled water.
4. Re-dyeing with Hematoxylin Solution for 5-8min and soak in tap water for 10min or redyeing with Methyl Green Solution for 2-3min.
5. View the sections under microscope after washing or sealing.

For Frozen Section

1. Restore the section to 37°C by immerse in water for 1-2min.
2. Dry in air and fix in TRAP Fixative for 1-3mins at 2-8°C.
3. Wash with distilled water and slightly dry in air.
4. Add the TRAP Incubation Solution to sections and place in 37°C incubator and incubate in dark for 45-60min. Then wash with distilled water.
5. Re-dyeing with Hematoxylin Solution for 3-5min and soak in tap water for 10min or redyeing with Methyl Green Solution for 2-3min.
6. View the sections under microscope after washing or sealing.

For Paraffin Section

1. Dewax paraffin sections by xylene twice for 5-10min each. Then rehydrate in graded ethanol (include 5min in 100% ethanol, 2min each in 90% ethanol and 70% ethanol).
2. Finally wash with distilled water for 2min, slightly dry in air.
3. Add the TRAP Incubation Solution to sections and place in 37°C incubator and incubate in dark for 45-60 min. Then wash with distilled water.

4. Re-dyeing with Hematoxylin Solution for 5-8min and soak in tap water for 10min or redyeing with Methyl Green Solution for 2-3min.
5. View the sections under microscope after washing or sealing.

Result

Positive	Purplish Red
Nucleus	Blue or Green

Clinical Significance

1. TRAP staining of hairy cells in hairy cell leukemia is strongly or moderately positive, and is not inhibited by tartaric acid. The other cells are negative or especially weak positive.
2. TRAP staining is positive in immature monocytes of acute leukemia, and weak positive in polymorphocytes. The response of progranulocytes to TRAP is different.
3. T lymphocyte TRAP staining is positive, with large and densely distributed granules. B lymphocyte is negative or weak positive with small granules.
4. Gaucher cells are strongly positive and Niemann-Pick cells are negative or weak positive.

Note

1. TRAP Incubation Solution is easy to lose effect. It's recommended to use skin puncture blood smears and stain quickly after drying.
2. When staining frozen sections, the exposure time of sections at room temperature should be reduced.
3. Reagent B2 is a storage solution for azo dyes, and the presence of granular suspensions is a normal phenomenon that does not affect dyeing. If you are concerned, pre mixing B1 and B2 can completely dissolve the particulate matter, but it needs to be used within 2 hours.
4. Samples should be fresh and treated immediately after sampling, otherwise the enzyme activity will be affected.
5. Tissue should be fixed in refrigerator at 2-8°C no more than 24 h, otherwise the activity of enzyme will be weakened or disappeared.
6. When embedding the tissue in paraffin, the temperature should not be higher than 56°C. It's recommended to use paraffin wax with melting point of 52-56°C for wax soaking. The soaking time should be short, otherwise the enzyme activity will weaken or disappear.
7. The impure dimethylbenzene will decompose the black precipitation, and the xylene above AR grade should be selected.
8. For your safety and health, please wear experimental clothes and disposable gloves.