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## Congo Red Amyloid Stain Kit(Highman Method)

**Cat:** G1535

**Size:** 3×50ml

**Storage:** RT, avoid light, valid for 6 months.

### Introduction

Amyloid is a kind of extracellular acidophilic substance with no fixed shape, which can exist in different tissues and organs, resulting in diseases called amyloidosis. Amyloid is mainly composed of proteins, most of which are arranged in reverse  $\beta$  - fold structure. Under the electron microscope, the amyloid materials are arranged as fibrils. In the case materials, there are a large number of non branching filaments, most of which are randomly arranged. The histological methods for the identification of amyloid substances include Violet Staining, Congo Red Staining and polarized light microscopy. In 1922, Bennhold found that Congo red can be used to identify starch like substances in vivo, and applied to tissue sections. After the improvement of highman, the dyeing effect is better.

### Kit Components

Reagent	3×50mL	Storage
Reagent (A): Highman Congo Red Staining Solution	50mL	RT, avoid light
Reagent (B): Alkaline Differentiation Solution	50mL	RT
Reagent (C): Lillie-Mayer Hematoxylin Staining Solution	50mL	RT, avoid light

### Self Provided Materials

10% formalin or distilled water, series of alcohol.

### Protocols(for reference only)

1. Conventionally fix in 10% neutral formalin, dehydrate and embed.
2. For paraffin section, dewax to distilled water; For frozen section, directly rinse into distilled water. If necessary, can remove the pigment.
3. Soak in Highman Congo Red Staining Solution and dye for 5-10mins.
4. Differentiate by Alkaline Differentiation Solution for 3-10s and then immediately remove to distilled water to stop differentiation. Wash with tap water.
5. Stain with Lillie-Mayer Hematoxylin Staining Solution for 1-2mins.
6. Slightly wash with tap water and replace the distilled water to make it differentiate and return blue.
7. Dehydrate by series of alcohol. Transparent by xylene. Seal with resinene.

### Result

Amyloid, Elastic Fiber, Eosinophilic Granules	Red
Nucleus	Blue

### Note

1. Section dewaxing should be as clean as possible, otherwise it will affect the dyeing effect.
2. The Alkaline Differentiation Solution should be kept in a closed state. Use it up as soon as possible once open.
3. The differentiation process of Alkaline Differentiation Solution is very important. It is necessary to enter water in time to stop differentiation to prevent excessive differentiation.
4. For your safety and health, please wear experimental clothes and disposable gloves.

### Related Products

G1534 Congo Red Amyloid Stain Kit(Modified Highman Method)

G1530 Congo Red Amyloid Stain Kit(Bennhold Method)

G1533 Congo Red Amyloid Stain Kit(Puchtler Method)

G1532 Congo Red Amyloid Stain Kit(Modified Stores Method)