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Collagen Fiber And Elastic Fiber Staining Kit(EVG-Verhoeff Method)

Cat: G1597

Size: 3×50mL/3×100mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent		3×50mL	3×100mL	Storage
Reagent(A) : Verhoeff Working Solution	Verhoeff Solution A	30mL	60mL	RT, avoid light
	Verhoeff Solution B	12mL	24mL	RT, avoid light
	Verhoeff Solution C	12mL	24mL	RT, avoid light
Before use, mix A,B and C as the ratio of 5:2:2 to prepare Verhoeff Working Solution. It is recommended to use in 2-3h.				
Reagent(B): Verhoeff Differentiation Solution		50mL	100mL	RT
Reagent(C): Modified VG Staining Solution		50mL	100mL	RT, avoid light

Introduction

Elastic fibers are mainly distributed in the arterial walls, alveolar walls, and skin of the human body. When fresh, they are yellow and have strong refractive properties. In routine tissue staining, it is often necessary to simultaneously display two types of fibers to study the coexistence relationship between fibers and their dependence on other structures. Based on this, an EVG composite staining kit was developed on the basis of a separate fiber staining kit, which is used in combination with conventional elastic fiber staining and VG collagen staining.

This reagent kit adopts a combined staining mechanism of Verhoeff Working Solution and modified VG staining solution, which can simultaneously stain elastic fibers and collagen fibers on a single slice without fading easily. After staining, the elastic fibers are purple black, the collagen fibers are red, and the muscles are yellow. It is suitable for observing whether there are abnormal lesions in the elastic fibers and collagen fibers in the tissue, thus assisting in diagnosis. This reagent kit has excellent quality, simple operation, stable staining, and beautiful and clear colors.

Protocol(for reference only)

1. Take fresh samples, fix and embed with paraffin. Cut the section into 3-8um.
2. After 50min of heat treatment on paraffin slices, dewax by xylene for 20min, then treat with gradient ethanol to distilled water for each step 5min.
3. Dye with Reagent(C): Modified VG Staining Solution for 10min, then wash with distilled water for 5-10s to remove excess dye.
4. Dye with Reagent(A): Verhoeff Working Solution for 5min. Wash with distilled water for 5-10s. (see note 1)
5. Differentiate in Reagent(B): Verhoeff Differentiation Solution for 5-10s, until the elastic fibers are clear. Wash with distilled water for 5-10s.
6. 75%, 85%, 95%, and 100% ethanol were subjected to conventional dehydration for 3-5 seconds each. Xylene was transparent twice for 1min each time, and neutral gum was used for sealing.
7. Microscopic observation of staining results.

Result

Elastic fiber	Purple Black
Collagen fiber	Red
Nucleus	Black

Note

1. The prepared reagent (A): Verhoeff Working Solution contains organic components, which are easy to diffuse and evaporate. During the staining process, it may be necessary to replenish the solution to prevent the slices from drying.
2. The degree of differentiation should be controlled under microscope to avoid over staining or light staining.
3. The Modified VG Staining Solution is stable in dyeing, not easy to fade after washing, and can be thoroughly cleaned.
4. For your safety and health, please wear laboratory clothes and disposable gloves.